

VENTFLO® ETCO₂/O₂ NASAL SAMPLING CANNULA

PROLONGED USE

LATEX FREE | SINGLE PATIENT USE | DISPOSABLE | DEHP FREE

The AirLife™ VentFLO® Extended EtCO₂/O₂ nasal sampling cannulas capture the exhaled CO₂ while simultaneously delivering oxygen. These EtCO₂/O₂ cannulas also feature a bright yellow reflective style connector, compatible with Microstream™, Capnostream™ and Oridion™ capnography monitors.



Microbial hydrophobic filter minimizes moisture ingress into monitor at the sampling port

3-Channel design throughout tubing and ear piece designed to help prevent kinking

Soft-Ears™ material helps ensure superior patient comfort, even during prolonged use

Nafion tubing allows the cannula to function longer in higher humidity conditions providing prolonged use


Yellow reflective style connector that indicates it is for extended use with ergonomic grip

Curved nares made of ultra-soft material

Fits-All O₂ connector comes standard on all sizes and lengths

VENTFLO® ETCO₂/O₂ NASAL SAMPLING CANNULA

- Compatible with capnography monitors that utilize Microstream™ technology, e.g. Medtronic Capnostream™ 35
- Simultaneously delivers oxygen while obtaining CO₂ sampling during spontaneous breathing
- Designed to prevent mixing of fresh oxygen with CO₂
- 3-Channel tubing prevents kinking
- Nafion tubing allows the cannula to function longer and in higher humidity conditions.
- Available in two tubing lengths, both with Fits-All O₂ connector
- Available as stand alone sampling line



ITEM	DESCRIPTION	LENGTH	PK
5107-SE-E	Adult Extended EtCO ₂ /O ₂ Nasal Cannula	7'	25
5114F-SE-E	Adult Extended EtCO ₂ /O ₂ Nasal Cannula	14'	25
5207F-SE-E	Pediatric Extended EtCO ₂ /O ₂ Nasal Cannula	7'	25


Microstream, Capnostream and Oridion are trademarks of Medtronic, related to products marketed by Medtronic.

VENTFLO® ETCO₂/O₂ ORAL-NASAL SAMPLING CANNULA

PROLONGED USE


LATEX FREE | SINGLE PATIENT USE | DISPOSABLE | DEHP FREE

The AirLife™ VentFLO® Extended EtCO₂/O₂ oral-nasal sampling cannulas capture the exhaled CO₂ while simultaneously delivering oxygen. These EtCO₂/O₂ cannulas also feature a bright yellow reflective style connector, compatible with Microstream™, Capnostream™ and Oridion™ capnography monitors.



VENTFLO® ETCO₂/O₂ ORAL-NASAL SAMPLING CANNULA

- Innovative design suitable for mouth-breathing patients
- Compatible with capnography monitors that utilize Microstream™ technology, e.g. Medtronic Capnostream™ 35
- Simultaneously delivers oxygen while obtaining CO₂ sampling during spontaneous breathing
- Designed to prevent mixing of fresh oxygen with CO₂
- Nafion tubing allows the cannula to function longer and in higher humidity conditions
- 3-Channel tubing prevents kinking
- Available in two tubing lengths, both with Fits-All O₂ connector



Microstream, Capnostream and Oridion are trademarks of Medtronic, related to products marketed by Medtronic.

ITEM	DESCRIPTION	LENGTH	PK
5707F-SE-E	Adult Extended EtCO ₂ /O ₂ Oral-Nasal Cannula	7'	25
5714F-SE-E	Adult Extended EtCO ₂ /O ₂ Oral-Nasal Cannula	14'	25

LIT-20051-C1 Rev. 1