


BINAPS® NASAL AIRFLOW PRESSURE AND SNORE TRANSDUCER

LATEX FREE | REUSABLE

BiNAPS® is a diagnostic pressure transducer device for use during sleep disorder studies. It will receive airflow/snore pressure signals via Salter Labs® Sleep Diagnostic Cannulas, process through a special piezo type pressure transducer and deliver an accurate signal to your polysomnograph input box. The exclusive integral lead verification tester has a green light to indicate functional leads from the BiNAPS to your input box.



BINAPS® NASAL AIRFLOW PRESSURE AND SNORE TRANSDUCER

- Detect airflow pressure and snore changes
- No battery required
- Compact and lightweight
- Integral circuit tester/indicator for lead verification
- Self adhering hook and loop mounting patch included

ITEM	DESCRIPTION	PK
5500	BiNAPS Nasal Airflow Pressure/Snore Transducer	1