



# TRACHEOSTOMY TUBE

LATEX FREE | STERILE | SINGLE USE | DISPOSABLE

AirLife® tracheostomy tubes provide features that aid in patient comfort from the low pressure cuff to the 360° swivel connector.



Obturator pre-loaded in tube



## TRACHEOSTOMY TUBE

- Obturator and tube clearly marked
- 15 mm Connector swivel connector
- Soft PVC
- Secure neck band retainer
- Cushion neck band comfortable for patient
- Radiopaque strip embedded for exact location of tube
- Low pressure cuff

## CUFFED

ITEM	SIZE	ID	LENGTH	PK
1-7393-50	20 Fr	5.0 mm	47 mm	10
1-7393-55	22 Fr	5.5 mm	50 mm	10
1-7393-60	24 Fr	6.0 mm	62 mm	10
1-7393-65	26 Fr	6.5 mm	67 mm	10
1-7393-70	28 Fr	7.0 mm	69 mm	10
1-7393-75	30 Fr	7.5 mm	75 mm	10
1-7393-80	32 Fr	8.0 mm	79 mm	10
1-7393-85	34 Fr	8.5 mm	83 mm	10
1-7393-90	36 Fr	9.0 mm	87 mm	10
1-7393-95	38 Fr	9.5 mm	91 mm	10
1-7393-10	40 Fr	10.0 mm	96 mm	10

## UNCUFFED

ITEM	SIZE	ID	LENGTH	PK
1-7390-30	12 Fr	3.0 mm	34 mm	10
1-7390-35	14 Fr	3.5 mm	37 mm	10
1-7390-40	16 Fr	4.0 mm	41 mm	10
1-7390-45	18 Fr	4.5 mm	43 mm	10
1-7390-50	20 Fr	5.0 mm	47 mm	10