

PARKER® SILICONE DISPOSABLE LARYNGEAL MASK

STERILE | LATEX FREE | SINGLE PATIENT USE | DISPOSABLE

The smooth contouring of our 100% medical grade silicone device enables easy insertion, is gentle to the airway, and is designed to produce an effective seal.



PARKER® DISPOSABLE LARYNGEAL MASK

- Smooth contouring of 100% silicone tube and cuff
- Non-fenestrated for easier insertion of diagnostic tools
- Not manufactured with natural rubber latex or DEHP
- Single patient use

ITEM	DESCRIPTION	SIZE	PK
H-PLMD-10-5	Disposable Laryngeal Mask	1	5
H-PLMD-15-5	Disposable Laryngeal Mask	1.5	5
H-PLMD-20-5	Disposable Laryngeal Mask	2	5
H-PLMD-25-5	Disposable Laryngeal Mask	2.5	5
H-PLMD-30-5	Disposable Laryngeal Mask	3	5
H-PLMD-40-5	Disposable Laryngeal Mask	4	5
H-PLMD-50-5	Disposable Laryngeal Mask	5	5
H-PLMDR-10-5	Wire Reinforced Disposable Laryngeal Mask	1	5
H-PLMDR-15-5	Wire Reinforced Disposable Laryngeal Mask	1.5	5
H-PLMDR-20-5	Wire Reinforced Disposable Laryngeal Mask	2	5
H-PLMDR-25-5	Wire Reinforced Disposable Laryngeal Mask	2.5	5
H-PLMDR-30-5	Wire Reinforced Disposable Laryngeal Mask	3	5
H-PLMDR-40-5	Wire Reinforced Disposable Laryngeal Mask	4	5
H-PLMDR-50-5	Wire Reinforced Disposable Laryngeal Mask	5	5

LIT-20061-C1 Rev. 1