

## VIXONE™ SMALL VOLUME JET NEBULIZER

LATEX FREE | SINGLE PATIENT USE | DISPOSABLE

The VixOne™ is a high-output nebulizer that delivers medication quickly and efficiently. Its versatile design and easy-to-fill cup with simple assembly makes this system a preferred solution for all patients.



VixOne™ Small Volume Jet Nebulizer

### VIXONE™ SMALL VOLUME JET NEBULIZER

- Popular acorn shape helps decrease medication waste by redirecting medication back to the jet
- Anti-spill “Tee” design helps keep exhaled condensation from contaminating medication
- Highly-efficient respirable mass output with MMAD of 2.7 microns promotes greater deposition
- Easy to fill cup with ¼ turn lock provides leak-proof seal
- Swivel mask feature allows treatment in upright positions or at angles up to 45 degrees
- 10 cc medication bowl capacity
- Very low dead volume (0.7 cc)
- Versatile design can be used in multiple configurations, including mouthpiece, mask, pacifier, reusable nebulizer cup, or in line with the ventilator



VixOne™ Small Volume Jet Nebulizer with Adult Mask



VixOne™ Small Volume Jet Nebulizer with Anti-spill “Tee”



PediNeb™ Nebulizer

# VIXONE™ SMALL VOLUME JET NEBULIZER

## VIXONE™ DISPOSABLE NEBULIZER

ITEM	DESCRIPTION	PK
0214	VixOne Disposable Nebulizer	50
0209	VixOne Disposable Nebulizer, Mouthpiece, Anti-spill "Tee", 7' Kink Resistant Tubing with Threaded Nut & Flexhose	50
0210	VixOne Disposable Nebulizer, Mouthpiece, Anti-spill "Tee", 7' Kink Resistant Tubing & Flexhose	50
0222	VixOne Disposable Nebulizer, Mouthpiece, Anti-spill "Tee", 7' Kink Resistant Tubing	50
0217	VixOne Disposable Nebulizer, Anti-spill "Tee", 7' Kink Resistant Tubing	50
0310	VixOne Disposable Nebulizer, with 7' Kink Resistant Tubing and Adult Mask	50
0311	VixOne Disposable Nebulizer, with 7' Kink Resistant Tubing and Pediatric Mask	50
0312	VixOne Disposable Nebulizer, with 7' Kink Resistant Tubing and Super Spike Pediatric Mask	50
0313	VixOne Disposable Nebulizer, with 7' Tubing, Threaded Nut and Super Spike Pediatric Mask	50
0319	VixOne Disposable Nebulizer, with 7' Tubing, Threaded Nut and Pediatric Mask	50

## VIXONE™ REUSABLE NEBULIZER

ITEM	DESCRIPTION	PK
0205	VixOne Reusable Nebulizer	10
0254	VixOne Reusable Nebulizer, Mouthpiece, Anti-spill "Tee", 7' Kink Resistant Tubing	10

## VIXONE™ NEBULIZER KITS

ITEM	DESCRIPTION	PK
0200	VixOne Multi Nebulizer Kit – Full Nebulizer Kit, 5 Extra Nebulizers	25

## PEDINEB™ NEBULIZER

ITEM	DESCRIPTION	PK
0382	VixOne Nebulizer, Pacifier, 45° Elbow & 7' Kink Resistant Tubing	12
0383	VixOne Nebulizer, Pacifier, Pediatric "Tee", Poppel Extension Set & 7' Kink Resistant Tubing	12
0385	Pedi-Neb, Pacifier with 45° Elbow	25