

INFU-SURG® PRESSURE INFUSION BAG

LATEX FREE | SINGLE PATIENT USE | DISPOSABLE | DEHP FREE | MR CONDITIONAL

The INFU-SURG® Pressure Infusion Bag is a safe, high-quality, reliable device for intravenous and intra-arterial infusion treatment, including A-line pressure monitoring. The Infu-Surg is a cost-effective, trusted product that medical professionals have relied on for decades. The dedicated fluid bag hook enables efficient loading and unloading of the system without having to remove the infusion bag from the IV pole.



360° view of the color-coded gauge facilitates accurate pressure monitoring while the safety pressure relief valve prevents over-inflation



Three-way stopcock ensures precise control of pressure while the large oval bulb allows for quick bladder inflation, all with ONE hand



Solution bag hook enables easy bag changing without having to remove the Infu-Surg from the IV pole; keeps fluid bag in proper position while in use

INFU-SURG® PRESSURE INFUSION BAG

- Pressure relief valve prevents over inflation (375 mm Hg pressure relief)
- Large, oval shaped bulb allows for quick and easy inflation of bladder
- Single-handed inflation and deflation design makes it easy to use and requires minimal training
- Dedicated fluid bag hook enables efficient loading and unloading of the system without having to remove the infusion bag from the IV pole
- Color-coded gauge makes for accurate pressure monitoring (0-300 mm Hg)
- Three-way stopcock ensures precise control of pressure
- Incredibly reliable – 100% tested

ITEM	DESCRIPTION	VOLUME	PK
4005H	Infu-Surg with Hook	500 mL	25
4010H	Infu-Surg with Hook	1000 mL	25
4030	Infu-Surg	3000 mL	25

LIT-20092-C1 Rev.1