


SALTER-STYLE® HIGH FLOW CANNULA

LATEX FREE | SINGLE PATIENT USE | DISPOSABLE

The comfort and convenience of oxygen therapy with nasal cannulas are now readily available for patients requiring higher oxygen flow rates. AirLife™ presents the High Flow Cannula from Salter Labs® that can replace cumbersome face masks with its enhanced comfort facepiece and 3-Channel safety tubing for flow rates up to 15 L/min. Extensive testing provides ample evidence of comparative performance.^{1,2,3}



FEATURES

- Delivers O₂ flows up to 15 L/min with 3-Channel safety tubing
- Patented slide bolo
- Enhanced reservoir facepiece
- Available in kits combining a high flow humidifier with a Salter Labs soft high flow cannula

PROBLEMS SOLVED

- Eliminates the need to use face masks
- Increased patient mobility
- Less cumbersome and claustrophobic than a mask
- Simple, secure adjustments

SALTER-STYLE® HIGH FLOW CANNULA

LATEX FREE | SINGLE PATIENT USE | DISPOSABLE

High Flow 1600HF Cannulas are anatomically designed and shaped with soft, curved nasal prongs to provide long-term comfort. Use the special 7900 High Flow Bubble Humidifier developed to partner with the High Flow 1600HF Cannulas, or order a High Flow Cannula Kit.

ADULT HIGH FLOW CANNULAS (FOR OXYGEN FLOWS UP TO 15 L/MIN)

| ITEM | DESCRIPTION | TUBING LENGTH | CONNECTOR | PK |
|-----------------|--|---------------|-----------------|----|
| 1600HF-7-25 | Salter-Style High Flow Cannula with Supply Tubing | 7' | Standard Ribbed | 25 |
| 1600HF-14-25 | Salter-Style High Flow Cannula with Supply Tubing | 14' | Standard Ribbed | 25 |
| 1600HF-25-10 | Salter-Style High Flow Cannula with Supply Tubing | 25' | Standard Ribbed | 10 |
| 16SOFT-HF-7-25 | Salter-Style Soft High Flow Cannula with Soft Tubing | 7' | Standard Ribbed | 25 |
| 16SOFT-HF-14-25 | Salter-Style Soft High Flow Cannula with Soft Tubing | 14' | Standard Ribbed | 25 |
| 16SOFT-HF-25-10 | Salter-Style Soft High Flow Cannula with Soft Tubing | 25' | Standard Ribbed | 10 |

ADULT HIGH FLOW CANNULA KITS (WITH SALTER-STYLE STANDARD OR SOFT DESIGN)

| ITEM | DESCRIPTION | TUBING LENGTH | CONNECTOR | PK |
|--------------|---|---------------|-----------------|----|
| 7907-7-10 | One (1) 1600HF-7 High Flow Cannula and One (1) 7900 Humidifier with Quad Nut | 7' | Standard Ribbed | 10 |
| 7914-14-10 | One (1) 1600HF-14 High Flow Cannula and One (1) 7900 Humidifier with Quad Nut | 14' | Standard Ribbed | 10 |
| 7925-25-10 | One (1) 1600HF-25 High Flow Cannula and One (1) 7900 Humidifier with Quad Nut | 25' | Standard Ribbed | 10 |
| 79SOFT-7-25 | One (1) 16SOFT-HF-7 High Flow Soft Cannula and One (1) 7900 Humidifier with Quad Nut | 7' | Standard Ribbed | 25 |
| 79SOFT-14-25 | One (1) 16SOFT-HF-14 High Flow Soft Cannula and One (1) 7900 Humidifier with Quad Nut | 14' | Standard Ribbed | 25 |
| 79SOFT-25-10 | One (1) 16SOFT-HF-25 High Flow Soft Cannula and One (1) 7900 Humidifier with Quad Nut | 25' | Standard Ribbed | 10 |

BUBBLE HUMIDIFIER

| ITEM | DESCRIPTION | CASE |
|-----------|--|------|
| 7900-0-25 | High Flow 350 cc Dry Bubble Humidifier for Flows From 6-15 L/min | 25 |



7900 Bubble Humidifier is simple to maintain and operate. Designed specifically for higher O₂ delivery.