

WATERPAP™ PRO

LATEX FREE | SINGLE PATIENT USE | DISPOSABLE | DEHP FREE

The WaterPAP™ PRO is a non-invasive device to deliver BCPAP (Bubble Continuous Positive Airway Pressure). BCPAP is a modality of respiratory support for preterm neonates with respiratory distress. CPAP helps splint the upper airways open, increasing the functional residual capacity while helping to decrease the work of breathing.



WATERPAP™ PRO

- Coaxial canister helps maintain a constant water level which improves accuracy with heated humidification
- Secure top provides protection to ensure tube stays in place
- Depth markers let health care providers quickly identify the CPAP setting
- Easy to change CPAP levels by adjusting tubing depth in water

ITEM	DESCRIPTION	CASE
WP-7800	WaterPAP Pro Setup – Coaxial Canister with 0-9 cm H ₂ O Graduated Lable, 10 mm Corrugated Tube Affixed with Air Diffuser, Tube Lock Top and Protective Cover	40
WP-7704	Bracket with Tube Guide	1