

CAPETTE™ AND STRAPETTE™ TUBE SUPPORTING HEADGEAR

LATEX FREE | SINGLE PATIENT USE | DISPOSABLE | DEHP FREE

The Capette™ is designed for the smallest neonate. The two-ply cotton hat helps to insulate the head and prevent heat loss. The Strapette™ and Tube Tackles are designed for use with neonatal and pediatric patients to help secure tubes without the use of adhesives. Together, Capette and Strapette provide a more effective tube support with less risk for skin injury associated with nasal cannula and nCPAP therapy.

- Compatible with nCPAP, high flow and low flow oxygen cannulas, and feeding tubes
- Adhesive free protecting neonates from skin injury
- Easy to apply soft foam head band that is easy to apply
- Tube tackles can attach anywhere to hold different types of tubing
- Head band and tube tackles can be easily adjusted in a matter of seconds for a more efficient assembly



Strapette with Small Tube Tackles



Medium Capette with Small Tube Tackles



Large Capette with Large Tube Tackles

CAPETTE™ TUBE SUPPORT

ITEM	DESCRIPTION	PATIENT	CASE
CP-4201	Small Capette with Small Tube Tackles	Neonates <800 g	40
CP-4202	Medium Capette with Small Tube Tackles	Neonates 800 - 1500 g	40
CP-4203	Large Capette with Small Tube Tackles	Neonates >1500 g	40
CP-4401	Small Capette with Large Tube Tackles	Neonates <800 g	40
CP-4402	Medium Capette with Large Tube Tackles	Neonates 800 - 1500 g	40
CP-4403	Large Capette with Large Tube Tackles	Neonates >1500 g	40



Small Tube Tackles



Large Tube Tackles

STRAPETTE™ TUBE SUPPORT

ITEM	DESCRIPTION	PATIENT	CASE
CP-4204	Strapette with Small Tube Tackles	Neonatal	40
CP-4220	Strapette with Small Tube Tackles	Pediatric	40
CP-4404	Strapette with Large Tube Tackles	Neonatal	40
CP-4420	Strapette with Large Tube Tackles	Pediatric	40

STRAPETTE™ TUBE SUPPORT

ITEM	DESCRIPTION	CASE
CP-0200	Small Tube Tackles, 3" (7 cm)	200
CP-0400	Large Tube Tackles, 4.5" (12 cm) for 10 mm CPAP Tubing	200