

SURGI-KIT®

LATEX FREE | SINGLE USE | DISPOSABLE

Surgi-Kit® Instrument Holders are disposable pockets for sterile delivery, common in surgical applications. This disposable tool maintains the sterile operating field and its simplicity eliminates cleaning procedures. Available sterile and non-sterile, choose from a variety of sizes and pocket configurations.


SURGI-KIT® INSTRUMENT AND SUTURE HOLDER

- Keeps instruments close to operating field
- Reduces time required for locating instruments
- Multipocket sleeves to accommodate a variety of instrument sizes
- Available with tape to adhere to operating drape
- Available in sterile and non-sterile options




Keep instruments within arm's reach by attaching instrument pocket to operating drape.

Select sizes available with peel-n-stick tape.




8 Pocket




4 Pocket




2 Pocket



8 Pocket + Tape



4 Pocket + Tape



2 Pocket + Tape



10 Pocket + Tape

STERILE

ITEM	DESCRIPTION	SIZE	PK
SK-100	8 pockets	19 cm x 23 cm	50
SK-150	4 pockets	19 cm x 11 cm	100
SK-200	2 pockets	28 cm x 18 cm	50
SK-300	8 pockets + tape	19 cm x 23 cm	50
SK-350	4 pockets + tape	19 cm x 11 cm	75
SK-400	2 pockets + tape	28 cm x 18 cm	50
SK-650	10 pockets + tape	5 cm x 19 cm	50

NON-STERILE

ITEM	DESCRIPTION	SIZE	PK
SK-100N	8 Pockets	19 cm x 23 cm	50
SK-150N	4 Pockets	19 cm x 11 cm	100
SK-200N	2 Pockets	28 cm x 18 cm	50
SK-300N	8 Pockets + Tape	19 cm x 23 cm	50
SK-350N	4 Pockets + Tape	19 cm x 11 cm	75
SK-400N	2 Pockets + Tape	28 cm x 18 cm	50
SK-650N	10 Pockets + Tape	5 cm x 19 cm	50

LAP-KIT®

LATEX FREE | SINGLE USE | DISPOSABLE

Lap-Kit® Laparoscopic Instrument Holders are disposable pockets for sterile delivery, common in surgical applications. Available as sterile or non-sterile, choose from a variety of sizes and pocket configurations.



LAP-KIT® LAPAROSCOPIC INSTRUMENT HOLDER


- Keeps laparoscopic instruments close to operating field
- Reduces time required for locating instruments
- Multipocket sleeves to accommodate a variety of instrument sizes
- Includes tape to adhere to operating drape
- Available in sterile and non-sterile options

STERILE


ITEM	DESCRIPTION	SIZE	PK
LK-700	2 Pockets + Tape	6.5 cm x 30.5 cm	25
LK-800	4 Pockets + Tape	3.2 cm x 30.5 cm	25
LK-900	3 Pockets + Tape	3.2 cm x 30.5 cm (2), 6.5 cm x 30.5 cm	25

NON-STERILE


ITEM	DESCRIPTION	SIZE	PK
LK-700N	2 Pockets + Tape	6.5 cm x 30.5 cm	25
LK-900N	3 Pockets + Tape	3.2 cm x 30.5 cm (2), 6.5 cm x 30.5 cm	25



2 Pockets



4 Pockets



3 Pockets