


MULTI-LINK™ X2 ECG LEADWIRES

LATEX FREE | SINGLE PATIENT USE

Multi-Link™ X2 is the simple, versatile, smart solution for ECG cables and leadwires. It allows standardization of ECG leadwires for all major monitoring platforms: GE™, Philips™, Mindray™, Nihon Kohden®, Drager™, Spacelabs™, Physio-Control™ and Zoll®.



Phillips set up with reusable 5-Lead set

GE set up with a disposable 5-lead set

MULTI-LINK™ ECG LEADWIRES

- Official, validated ECG supplies and accessories for GE Healthcare® patient monitoring systems
- Compatible with major bedside monitoring platforms including GE HealthCare, Drager, Mindray, Nihon Kohden and Spacelabs
- One patient, one leadwire
- Small and lightweight design maximizes patient experience
- May help reduce cross-contamination and associated hospital-acquired infections (HAI)
- May help patient workflow and room turnover time
- Tested for chemical resistance of cleaning and disinfection agents
- Compatibility tested within the hospital environment for signal accuracy according to AAMI 60601-2-27
- Leadwire lengths make it possible to utilize the same leadwires across different departments such as ER, OR, ICU and telemetry
- Defibrillation proof

ITEM	DESCRIPTION	SIZE	PK
2052133-005	Leadwire, ECG, 3-Lead, Grabbers, AHA	100cm/39"	10
2052133-027	Leadwire, ECG, 3-Lead, Grabbers, AHA	75cm/29"	10
2052133-006	Leadwire, ECG, 5-Lead, Grabbers, AHA	100cm/39"	10
2052133-007	Leadwire, ECG, 5-Lead, Grabbers, AHA	75cm/29"	10
2086367-003	Leadwire, ECG, 6-Lead, Grabbers, AHA	100cm/39"	10
2086367-004	Leadwire, ECG, 6-Lead, Grabbers, AHA	75cm/29"	10
2052133-008	Leadwire, ECG, V-Lead, Grabbers, AHA	100cm/39"	10

Not an all inclusive list of products. Contact your AirLife representative for a full list of offerings including snap connections and 12-ft length leadwires.