

SUPERNO₂VA™ ET MASK AND SYSTEM


LATEX FREE | SINGLE PATIENT USE | DISPOSABLE | DEHP FREE | MRI COMPATIBLE

The SuperNO₂VA™ Et Mask is a non-invasive positive pressure system designed to deliver oxygen, air and anesthetic gas to patients who have been diagnosed with or are at risk of obstructive sleep apnea.



Adult Medium Mask (SEM-20)


Adult Large Mask (SEL-20)



PRE-OP
Used to maintain O₂ saturation



INTRA-OP
Deep sedation and intra-oral



POST-OP
Provides positive pressure or passive oxygenation

SUPERNO₂VA ET

- Allows for effective pre-oxygenation and rescue ventilation without the need for additional equipment
- Maintains oral access
- Can provide passive oxygenations, high flow oxygen, and oral and/or nasal CO₂ capture.

ITEM	DESCRIPTION	SIZE	PATIENT HEIGHT	PK
SEM-20	SuperNO ₂ VA™ Et Mask	Adult Medium	4'7" - 6'9"	20
SEL-20	SuperNO ₂ VA™ Et Mask	Adult Large	5'10" - 7'3"	20
SSEM-20	SuperNO ₂ VA™ Et System	Adult Medium	4'7" - 6'9"	20
SSEL-20	SuperNO ₂ VA™ Et System	Adult Large	5'10" - 7'3"	20