


NONINVASIVE VENTILATION MASKS

LATEX FREE | DEHP FREE

AirLife™ noninvasive ventilation (NIV) full face and nasal masks offer clinicians a reliable option while providing patients comfort and security. The adjustable four-point headgear enables ease in fit and adjustment while the flexible swivel elbow reduces torque applied to the seal between the mask and patient, helping to minimize pressure points.



Adjustable forehead rest for comfort and stability

Flexible swivel connector enables ease in positioning of the circuit



Quick-release clip eliminates the need to refit after removal



Nasal Mask

MASK FLEXES WHEN THE CIRCUIT IS PULLED



Mask when flexed



Mask not flexed

FEATURES

- Comfortable four-point headgear with crown strap for better stability and sealing
- Soft seal that conforms to the face, reducing need to overtighten
- Lightweight, silicone design that contours to the patient's face for a comfortable fit
- Include both vented and non-vented options and are available in three sizes to fit a wide range of patients.

AIRLIFE™ NIV MASKS

ITEM	DESCRIPTION	SIZE	PK
NIV040S	Full Face Mask, Vented with Anti-Asphyxia Valve (AAV)	Small	10
NIV040M	Full Face Mask, Vented with AAV	Medium	10
NIV040L	Full Face Mask, Vented with AAV	Large	10
NIV043S	Full Face Mask, Non-Vented with AAV	Small	10
NIV043M	Full Face Mask, Non-Vented with AAV	Medium	10
NIV043L	Full Face Mask, Non-Vented with AAV	Large	10
NIV041S	Full Face Mask, Non-Vented	Small	10
NIV041M	Full Face Mask, Non-Vented	Medium	10
NIV041L	Full Face Mask, Non-Vented	Large	10
NIV042S	Nasal Mask, Vented	Small	10
NIV042M	Nasal Mask, Vented	Medium	10
NIV042L	Nasal Mask, Vented	Large	10
NIV044S	Nasal Mask, Non-Vented with AAV	Small	10
NIV044M	Nasal Mask, Non-Vented with AAV	Medium	10
NIV044L	Nasal Mask, Non-Vented with AAV	Large	10

LIT-20113-C1 Rev1