

TEMPERATURE MONITORING

LATEX FREE | REUSABLE

AirLife™ reusable temperature cables are validated accessories for use with GE HealthCare® bedside patient monitors. These durable reusable cables can be utilized with both disposable and reusable temperature probes



2016998-001



2021700-001



165640



2021701-001



2021704-002

TEMPERATURE CABLES

- For use with GE HealthCare patient monitors
- Connects directly to reusable probes or to interconnect cable (165640 and 165641) for disposable probes

ITEM	DESCRIPTION	CABLE LENGTH	PK
2016998-001	Dual Temperature Cable	0.5 m (1.7')	1
2021700-001	Single Temperature Cable	0.5 m (1.7')	1

INTERCONNECT CABLE

- Interconnect cable for disposable probes
- For use with 2016998-001 or 2021700-001

ITEM	DESCRIPTION	CABLE LENGTH	PK
165640	Interconnect Cable	1.3 m (4.3')	1
165641	Interconnect Cable	2.8 m (9')	1

DIRECT CONNECT CABLE

- Connects disposable probes directly to GE HealthCare patient monitors
- No interconnect cable required

ITEM	DESCRIPTION	CABLE LENGTH	PK
2021701-001	Direct Connect Cable, Single	3.6 m (12')	1
2021704-002	Direct Connect Cable, Dual	3.6 m (12')	1