

# STANDARD ANESTHESIA FACE MASK

LATEX FREE | DISPOSABLE

Standard disposable face masks - featuring color-coded hook rings, brilliant clarity and semi-soft teardrop face plates.



## FEATURES

- A dual inflation valve (tail valve) permits cushion adjustment with a syringe or by mouth.
- Soft, pliable, air-filled cushion molds easily to face to deliver a low-pressure seal.

## STANDARD ANESTHESIA FACE MASKS

ITEM	DESCRIPTION	SIZE	PK
4810VS	Standard, Top Valve, Air Filled Cushion Mask, Red Ring	1 Neonate	20
4820VS	Standard, Top Valve, Air Filled Cushion Mask, Green Ring	2 Infant	20
4830VS	Standard, Top Valve, Air Filled Cushion Mask, Purple Ring	3 Toddler	20
4840VS	Standard, Top Valve, Air Filled Cushion Mask, Yellow Ring	4 Child	20
4852VS	Standard, Top Valve, Air Filled Cushion Mask, White Ring	5 Small Adult	100
4860VS	Standard, Top Valve, Air Filled Cushion Mask, Blue Ring	6 Medium Adult	20
4861VS	Standard, Top Valve, Air Filled Cushion Mask, Blue Ring	6 Medium Adult	100
4870VS	Standard, Top Valve, Air Filled Cushion Mask, Orange Ring	7 Large Adult	20

LIT-20119-C1 Rev1