

ROUND ANESTHESIA FACE MASK

LATEX FREE | DISPOSABLE

These traditional line of non-flexible, rounded dome face masks provides several styles to suit individual needs.



FEATURES

- A dual inflation valve (tail valve) permits cushion adjustment with a syringe or by mouth.
- Soft, pliable, air-filled cushion molds easily to face to deliver a low-pressure seal.

ROUND ANESTHESIA FACE MASKS

ITEM	DESCRIPTION	SIZE	PK
5210	Round Rigid Dome, Air Filled Cushion, Clear w/o Valve	1 Neonate	20
5211	Round Rigid Dome, Air Filled Cushion, Clear, Tail Valve	1 Neonate	100
5220	Round Rigid Dome, Air Filled Cushion, Clear w/o Valve	2 Infant	20
5221EU	Round Rigid Dome, Air Filled Cushion, Clear w/o Valve	2 Infant	100
5230	Round Rigid Dome, Air Filled Cushion, Clear, Tail Valve	3 Toddler	20
5231EU	Round Rigid Dome, Air Filled Cushion, Clear, Tail Valve	3 Toddler	100
5240	Round Rigid Dome, Air Filled Cushion, Clear, Tail Valve	4 Child	20
5241	Round Rigid Dome, Air Filled Cushion, Clear, Tail Valve	4 Child	100
5250EU	Round Rigid Dome, Air Filled Cushion, Clear, Tail Valve	5 Adult	20
5245	Rigid Dome, Air Filled Cushion, Clear, Tail Valve	5 Adult	20
5246	Rigid Dome, Air Filled Cushion, Clear, Tail Valve	5 Adult	100
5251VS	Round Rigid Dome, Air Filled Cushion, Clear, Tail Valve	5 Adult	100

LIT-20149-C1 Rev2