


UNI-LIM™ AND PEDI-LIM™ ANESTHESIA CIRCUIT

LATEX FREE | SINGLE PATIENT USE | DISPOSABLE


The single-limb circuit is the optimal choice for your anesthesia circuit needs. It is designed to function as a standard circuit anesthesia breathing system with streamlined tube-within-a-tube and eliminated bulky connectors for a simple, efficient anesthesia circuit.




HEALTHTRUST®
Contract #102574
HCA Formulary



9810 Anesthesia Circuit Kit, 72" Uni-lim coaxial circuit, patient elbow with capped gas sample port, M/M GSL 10' , 2 B/V filter, 3L latex-free breathing bag, disposable




9812 Anesthesia Circuit Kit, 60" Pedi-lim coaxial circuit, patient elbow with capped gas sample port, M/M GSL 10' , 2 B/V filter, 1L latex-free breathing bag, disposable



UNI-LIM™ AND PEDI-LIM™ ANESTHESIA CIRCUITS

- Features a streamlined tube-within-a-tube design which eliminates twisting between the inspiratory and expiratory limbs
- Eliminates bulky connectors at the patient interface
- Provides natural heat and moisture conservation properties to that of a standard circle circuit
- Engineered with a secure-t inspiratory limb to prevent patient end disconnects — no floating internal tubes
- ISO-standard connections for secure, leak resistant attachments



UNI-LIM™ AND PEDI-LIM™ ANESTHESIA CIRCUIT

LATEX FREE | SINGLE PATIENT USE | DISPOSABLE

CURRENT HEALTHTRUST (AMBU) SKU	AIRLIFE CROSS	DESCRIPTION	UOM	UOM FACTOR
6753F-19021Z	9239	Anesthesia Circuit Kit, 72" expandable dual limb circuit, patient elbow with capped gas sample port, Size 6 large adult mask, M/M GSL 10', 2 B/V filter, 3L latex-free breathing bag, disposable	CA	20
8903F-61Z	A5Z12XXX	Anesthesia Circuit Kit, 96" expandable dual limb circuit, patient elbow with capped gas sample port, 1 B/V filter, 3L latex-free breathing bag, disposable	CA	20
D365-6121Z	9803	Anesthesia Circuit Kit, 60" Uni-lim coaxial circuit, patient elbow with capped gas sample port, Size 5 adult mask, M/M GSL, 2 B/V filter, 3L latex-free breathing bag, disposable	CA	20
D366-6121Z	9804	Anesthesia Circuit Kit, 60" Uni-lim coaxial circuit, patient elbow with capped gas sample port, Size 6 large adult mask, M/M GSL 10', 2 B/V filter, 3L latex-free breathing bag, disposable	CA	20
D460-8021Z	9169	Anesthesia Circuit Kit, 60" Uni-lim coaxial circuit, patient elbow with capped gas sample port, M/M GSL 10', 2 B/V filter, 2L latex-free breathing bag, disposable	CA	20
DF3110-6121Z	9802	Anesthesia Circuit Kit, 120" Uni-lim coaxial circuit, patient elbow with capped gas sample port, M/M GSL 10', 2 B/V filter, 3L latex-free breathing bag, disposable	CA	20
DF3115-6121Z	9806	Anesthesia Circuit Kit, 120" Uni-lim coaxial circuit, patient elbow with capped gas sample port, Size 5 adult mask, M/M GSL, 2 B/V filter, 3L latex-free breathing bag, disposable	CA	20
DF370-6121Z	9810	Anesthesia Circuit Kit, 72" Uni-lim coaxial circuit, patient elbow with capped gas sample port, M/M GSL 10', 2 B/V filter, 3L latex-free breathing bag, disposable	CA	20
DF370-61Z	9298	Anesthesia Circuit Kit, 72" Uni-lim coaxial circuit, patient elbow with capped gas sample port, 2 B/V filter, 3L latex-free breathing bag, disposable	CA	20
DF375-17621Z	9809	Anesthesia Circuit Kit, 72" Uni-lim coaxial circuit, patient elbow with capped gas sample port, sampling tee, Size 5 adult mask, M/M GSL 10', 2 B/V filter, 3L latex-free breathing bag, disposable	CA	20
DF375-6121Z	9694	Anesthesia Circuit Kit, 72" Uni-lim coaxial circuit, patient elbow with capped gas sample port, Size 5 adult mask, M/M GSL 10', 2 B/V filter, 3L latex-free breathing bag, disposable	CA	20
DF470-6121Z	9805	Anesthesia Circuit Kit, 72" Uni-lim coaxial circuit, patient elbow with capped gas sample port, M/M GSL 10', 2 B/V filter, 2L latex-free breathing bag, disposable	CA	20
DF475-6121Z	9177	Anesthesia Circuit Kit, 72" Uni-lim coaxial circuit, patient elbow with capped gas sample port, Size 5 adult mask, M/M GSL, 2 B/V filter, 2L latex-free breathing bag, disposable	CA	20
KN6903B-6121Z	9665	Anesthesia Circuit Kit, 96" expandable dual limb circuit, patient elbow with capped sampling port, in-line mini HMEF with sample port, Size 6 large adult mask, M/M GSL 10', 3L latex-free breathing bag, all disposable	CA	20
KNCF376-6121Z	9808	Anesthesia Circuit Kit, 72" Uni-lim coaxial circuit, patient elbow with capped gas sample port, Size 6 large adult mask, in-line mini HMEF with sample port, M/M GSL 10', 3L latex-free breathing bag, disposable	CA	20
KNUF3116-6121Z	9811	Anesthesia Circuit Kit, 120" Uni-lim coaxial circuit, patient elbow with capped gas sample port, in-line HMEF with sampling port, Size 6 large adult mask, M/M GSL 10', 3L latex-free breathing bag, disposable	CA	20
KNUF4116-6121Z	9807	Anesthesia Circuit Kit, 120" Uni-lim coaxial circuit, patient elbow with capped gas sample port, in-line HMEF with sampling port, Size 6 large adult mask, M/M GSL 10', 2L latex-free breathing bag, disposable	CA	20
PDP160-6121	9812	Anesthesia Circuit Kit, 60" PEDI-lim coaxial circuit, patient elbow with capped gas sample port, M/M GSL 10', 2 B/V filter, 1L latex-free breathing bag, disposable	CA	20

LIT-20167-C1 Rev. 1