

MX40 ADAPTER

LATEX FREE

The Multi-Link™ MX40 Adapter connects to the Philips™ IntelliVue MX40 telemetry device and is designed for use with AirLife's Multi-Link X2 leadwires.



FEATURES

- SpO₂ port features an effortlessly accessible opening with a securing mechanism preventing the SpO₂ device from disconnecting
- Compatible with Masimo LNCS pulse oximeter sensors and Philips 9-pin pulse oximeter sensors
- Available for both 5 lead and 6 lead
- Tested for chemical resistance of cleaning and disinfection agents instructed in product IFU

ITEM	DESCRIPTION
2090101-101	Tele Adapter X2, 5-Lead for Philips MX40*
2090101-102	Tele Adapter X2, 6-Lead for Philips MX40