



AIRLIFE OXYGEN CANNULAS

LATEX FREE | DISPOSABLE


The AirLife line of oxygen cannulas offer a premium selection of nasal cannulas designed for patient well-being and care. This line features soft head tubing to minimize chafing, clear cushion designs, and foam ear covers for reduced skin irritation. Additionally, our standard nasal cannulas include angulated lip plates and beveled tips for enhanced comfort.



Soft Head Tubing Nasal Cannula




Cushion Nasal Cannula



Foam Cover Ear Cushions Nasal Cannula

AIRLIFE OXYGEN CANNULAS

- Soft head tubing material for improved patient comfort
- Crush-resistant tubing to minimize occlusion and kinking
- Clear cushion design for optimal patient comfort
- Foam ear covers for reduced skin irritation
- Angulated, flexible lip plate with beveled tips for softer edge
- Available in sizes from neonate to adult with varying lengths of oxygen tubing
- Universal connector to eliminate the need for a christmas tree adapter*



Standard, Economical - Flared

Standard, Economical - Non-Flared

ITEM	DESCRIPTION	SIZE	LENGTH	PK
SFT2699	Soft Head Tubing	Adult	7'	50
SFT2614U	Soft Head Tubing	Adult	14'	50
SFT2625U	Soft Head Tubing	Adult	25'	25
SFT2692	Soft Head Tubing	Pediatric	7'	50
SFT2601	Soft Head Tubing	Infant	7'	50
SFT2611	Soft Head Tubing	Neonate	7'	50
002699	Cushion Face Piece	Adult	7'	50
002600-14U	Cushion Face Piece	Adult	14'	50
002600-25U	Cushion Face Piece	Adult	25'	25
002692	Cushion Face Piece	Pediatric	7'	50
002601	Cushion Face Piece	Infant	7'	50
002611	Cushion Face Piece	Neonate	7'	50
FM2699	Foam Cover Ear Cushions	Adult	7'	50
FM2614U	Foam Cover Ear Cushions	Adult	14'	50
FM2625U	Foam Cover Ear Cushions	Adult	25'	25
FM2692	Foam Cover Ear Cushions	Pediatric	7'	50
001365	Standard, Economical - Flared	Adult	7'	50
001321U	Standard, Economical - Flared	Adult	14'	50
001322U	Standard, Economical - Flared	Adult	21'	25
001366	Standard, Economical - Non-Flared	Adult	7'	50
001326U	Standard, Economical - Non-Flared	Adult	14'	50
001327U	Standard, Economical - Non-Flared	Adult	21'	25
001329U	Standard, Economical - Non-Flared	Adult	50'	15

* Universal connector not available on infant and neonate size cannulas