

# AirFlow™



Pictured Adult AF5140MB

## The SunMed® AirFlow™ Resuscitation Bag offers new, enhanced features for higher levels of accuracy and safety.

- Redesigned manometer to improve visibility and enhance accuracy.
- Improved accuracy on our patented, integrated color-coded pressure manometer.<sup>1</sup>
- MRI Conditional up to 7 tesla.
- Smaller patient valve allows for decreased dead space < 8 mL.<sup>2</sup>
- Sampling port allows for EtCO<sub>2</sub> monitoring.

<sup>1</sup>.  $\pm 3$  cm H<sub>2</sub>O Up to 15cm H<sub>2</sub>O  
<sup>2</sup>. < 5mL + 10% of delivered volume



Scan to watch  
AirFlow video