


COMFORT SOFT PLUS® OXYGEN CANNULAS

The unique elastomeric properties of the over-the-ear Comfort Soft Plus tubing minimizes friction and shear forces on the skin at all contact points and acts as a “strain relief” when a patient turns their head. As a result, the AirLife® Comfort Soft Plus cannulas are more comfortable, enhancing patient compliance and therapy outcomes.



COMFORT SOFT PLUS® OXYGEN CANNULAS

ITEM	DESCRIPTION	SIZE	TUBE LENGTH	CONNECTOR	CASE
0567	Comfort Soft Plus Oxygen Cannula	Adult	7'	Threaded Nut	50
0153	Comfort Soft Plus Oxygen Cannula	Premature	7'	Threaded Nut	25
0163	Comfort Soft Plus Oxygen Cannula	Neonatal	7'	Threaded Nut	25
0173	Comfort Soft Plus Oxygen Cannula	Infant	7'	Threaded Nut	50
0183	Comfort Soft Plus Oxygen Cannula	Pediatric	7'	Threaded Nut	50

PS-20026-C1 Rev1