

User Guide

8900 SERIES DISPOSABLE SMALL VOLUME NEBULIZER KIT

Intended Use

The 8900 Small Volume Nebulizer is a hand-held pneumatically powered nebulizer used to deliver aerosolized medication via a mouthpiece or an aerosol mask. The patient population includes infant, pediatric and adult patients. The nebulizer is a disposable, single-patient, multiple-use medical device.

Caution

Do not reuse on more than one patient.
Do not sterilize the device.

Contraindications

No known contraindications.

Packaged

The 8900 Small Volume Nebulizer kit is individually packaged with a:

1. Mouthpiece or aerosol mask
2. Tee adapter
3. 6 inch aerosol flex tube, a.k.a. reservoir tube
4. Nebulizer (cap, cone and medication cup)
5. 3-Channel supply tubing, 7 feet

Note: The 8900 components are available separately and packaged in different configurations.

