

**00**X Small
Neonate

Vent Mask II™

Non-Inflatable Non-Conductive Anesthesia Mask

REF VR0040 / 00 / X Small Neonate**Directions for Use**

Place mask over patient's nose and mouth to direct gases to upper airway.

**Caution**

Do not attempt to sterilize, autoclave or clean this device.



Single Use

**NON
STERILE****LATEX****PHT****DEHP****R_X** ONLYVentlab, LLC
2710 Northridge Dr. NW, Suite A
Grand Rapids, MI 49544 USAManufactured for Ventlab, LLC,
a Subsidiary of SunMed Holdings, LLC
www.Sun-Med.com
Made in China by SunMed AP, LTD

VR0040.IFU Rev.3 2022/11