



SUNMED® LARYNGEAL MASK AIRWAY

LATEX FREE | SINGLE PATIENT USE | DISPOSABLE | DEHP FREE | STERILE

The SunMed® Laryngeal Mask Airways are designed to help achieve and maintain airway control during emergency and routine procedure in fasted patients and difficult airway situations. Both Standard and Flex variations function as less-invasive alternatives to the classic endotracheal tube. These products provide a low-cost option while still delivering the known benefits of a supraglottic airway device.



Color coded tabs support quick & easy identification



SUNMED® LARYNGEAL MASK AIRWAYS

ITEM	AIRWAY TYPE	MASK SIZE	PATIENT WEIGHT	MAX CUFF INFLATION VOLUME	TAB COLOR
1-2021-01	Standard	1	<11 lbs (<5 kg)	4 mL	Pink
1-2021-15	Standard	1.5	11-22 lbs (5-10 kg)	7 mL	Blue
1-2021-02	Standard	2	22-44 lbs (10-20kg)	10 mL	Black
1-2021-25	Standard	2.5	44-66 lbs (20-30 kg)	14 mL	White
1-2021-03	Standard	3	66-110 lbs (30-50 kg)	20 mL	Green
1-2021-04	Standard	4	110 - 154 lbs (50-70 kg)	30 mL	Yellow
1-2021-05	Standard	5	154-220 lbs (70-100 kg)	40 mL	Red
1-2021-06	Standard	6	>220 lbs (>100 kg)	50 mL	Blue
1-0372-10	Flex	1	<11 lbs (<5 kg)	4 mL	Pink
1-0372-15	Flex	1.5	11-22 lbs (5-10 kg)	7 mL	Blue
1-0372-20	Flex	2	22-44 lbs (10-20kg)	10 mL	Black
1-0372-25	Flex	2.5	44-66 lbs (20-30 kg)	14 mL	White
1-0372-30	Flex	3	66-110 lbs (30-50 kg)	20 mL	Green
1-0372-40	Flex	4	110 - 154 lbs (50-70 kg)	30 mL	Yellow
1-0372-50	Flex	5	154-220 lbs (70-100 kg)	40 mL	Red
1-0372-60	Flex	6	>220 lbs (>100 kg)	50 mL	Blue

LIT-20200-C1 Rev.1