

AIRLIFE DUOTHERM® HUMIDIFICATION SYSTEM

SINGLE USE | NON-STERILE | LATEX FREE | DEHP FREE

The AirLife DuoTherm® humidification system is designed as an integral part of the personalized care environment. By providing improved clinical and workflow efficiencies, the AirLife DuoTherm System helps clinicians personalize humidification delivery to fit the needs of each patient over the broad spectrum of patient populations.

FEATURES

- Greater adjustability than other systems to personalize humidification delivery to all patient types (neonatal through adult)
- Delivery of therapeutic levels of humidification at greater gas flow ranges than other similar systems
- Compatible with invasive and noninvasive therapy systems commonly found within institutional care environments
- Clean and clear display of AirLife DuoTherm operating status and personalized humidification delivery including alerts to clinicians if conditions are outside of the desired operating range
- Dual-walled insulated circuits that remain stable even as environmental conditions change



HUMIDIFICATION SYSTEM & CONSUMABLES

ITEM	DESCRIPTION	PK
377HTR	AirLife DuoTherm Humidification Heater	1
377CBRE	AirLife DuoTherm Humidification Chamber	10
377PWR	AirLife DuoTherm 110V Power Adapter	1
377ASL	AirLife DuoTherm Heater Wire Adapter, Reusable, Single Limb	1
377ADL	AirLife DuoTherm Heater Wire Adapter, Reusable, Dual Limb	1
377PR5E	AirLife DuoTherm 5' Reusable Temperature Probe	1
377PR6E	AirLife DuoTherm 6' Reusable Temperature Probe	1
37712CE	AirLife DuoTherm Neonate Single Limb Circuit Kit, Dual-Wall Heated Wire Circuit with Chamber	10
37723CE	AirLife DuoTherm Neonate Dual Limb Circuit Kit, Dual-Wall Heated Wire Circuit with Chamber	10
37736CE	AirLife DuoTherm Pediatric Single Limb Circuit Kit, Dual-Wall Heated Wire Circuit with Chamber	10
37748CE	AirLife DuoTherm Pediatric Dual Limb Circuit Kit, Dual-Wall Heated Wire Circuit with Chamber	10
37756CE	AirLife DuoTherm Adult Single Limb NIV Circuit Kit, Large Bore Heated Wire Circuit with Chamber	10
37747CE	AirLife DuoTherm Adult Dual Limb Circuit Kit, Dual-Wall Heated Wire Circuit with Chamber	10
37735CE	AirLife DuoTherm Adult Single Limb Circuit Kit, Dual-Wall Heated Wire Circuit with Chamber	10
37740CE	AirLife DuoTherm Adult Single Limb NIV Circuit Kit, Dual-Wall Heated Wire Circuit with Chamber	1
AH115-05E	Pressure Relief Valve, Infant, 40 cm H ₂ O	5

COMPONENTS



BRACKETS

ITEM	DESCRIPTION	PK
377BR101	Bracket, Pole Clamp Mounting	1
377BR102	Bracket, Mounting	1
377BR103	Bracket, Pole & Rail Mounting	1
377BR201	Water Pole	1
377BR301	Bracket, Offset Pole & Rail Mounting	1
377BR302	Bracket, PB840/980 Mounting	1
377BR303	Bracket, 44 mm, Adapter & Mounting	1
377BR304	Bracket, Draeger, Vent-Mounting	1
377BR305	Bracket, Mounting, with Adapter – Concha	1