

BALLARD™ NEONATAL AND PEDIATRIC CLOSED SUCTION SYSTEMS

MR CONDITIONAL

The Ballard™ Closed Suction Systems for Neonates and Pediatrics are designed to safely suction patients on mechanical ventilation, while providing protection for the patient and the care-giver.



CLOSED SUCTION FEATURES

- Low durometer catheter and rounded atraumatic tip to minimize tracheal trauma
- Elbow and “Y” configurations for standard and high-frequency ventilation
- Available in 5, 6, 7, 8, 10 & 12 French sizes
- Reusable up to 24 hours
- WET PAK configurations include 3 mL saline bullets
- ET Tube Adapters available in sizes 2 – 4.5 mm

Available Systems: “Y”, Elbow, Elbow Manifold, and Multi-Access Catheters for surfactant delivery.

BALLARD™ NEONATAL AND PEDIATRIC CLOSED SUCTION SYSTEMS

MR CONDITIONAL

CLOSED SUCTION SYSTEMS -“Y”

ITEM	SIZE	LENGTH	ADAPTERS	PK
195	5 Fr	30.5 cm / 12"	Y-Endotracheal Tube Adapters 2 mm, 2.5 mm	20
196	6 Fr	30.5 cm / 12"	Y-Endotracheal Tube Adapters 2.5 mm, 3 mm, 3.5 mm	20
197	7 Fr	30.5 cm / 12"	Y-Endotracheal Tube Adapters 3.0 mm, 3.5 mm, 4.0 mm	20
198	8 Fr	30.5 cm / 12"	Y-Endotracheal Tube Adapters 3.0 mm, 3.5 mm, 4.0 mm	20
1910	10 Fr	40.5 cm / 16"	Y-Endotracheal Tube Adapters 4.5 mm, 5.0 mm, 5.5 mm	20
1912	12 Fr	40.5 cm / 16"	Y-Endotracheal Tube Adapters 5.5 mm, 6.0 mm, 6.5 mm	20



CLOSED SUCTION SYSTEM - “Y” WITH SALINE VIALS

ITEM	SIZE	LENGTH	ACCESSORIES	PK
1968127	6 Fr	30.5 cm / 12"	Y- Adapters (2.5 mm, 3 mm, 3.5 mm), 1 Oral Suction System, 1 Package Suction Probes, 1 9' Bubble Tubing, 6 3 mL Saline Vials	20
1978127	7 Fr	30.5 cm / 12"	Y- Adapters (3 mm, 3.5 mm, 4 mm), 5 8 Fr Suction Probes, 1 9' Bubble Tubing, 1 Parallel Y, 1 Connection tubing with suction valve, 6 3 mL Saline Vials	20
1988128	8 Fr	30.5 cm / 12"	Y- Adapters (3 mm, 3.5 mm, 4 mm), 1 Oral Suction System, 1 Package Suction Probes, 1 Parallel Y, 1 9' Bubble Tubing, 6 3mL Saline Vials	20



CLOSED SUCTION SYSTEMS - ELBOW

ITEM	SIZE	LENGTH	ADAPTERS	PK
202	5 Fr	30.5 cm / 12"	Endotracheal Tube Adapters 2 mm, 2.5 mm	20
206	6 Fr	30.5 cm / 12"	Endotracheal Tube Adapters 2.5 mm, 3 mm, 3.5 mm	20
207	7 Fr	30.5 cm / 12"	Endotracheal Tube Adapters 2.5 mm, 3 mm, 3.5 mm, 4 mm	20
208	8 Fr	35.5 cm / 14"	Endotracheal Tube Adapters 3 mm, 3.5 mm, 4 mm	20
210	10 Fr	40.5 cm / 16"	—	20



CLOSED SUCTION SYSTEMS - ELBOW MANIFOLD

ITEM	SIZE	LENGTH	ADAPTERS	PK
200	6 Fr	30.5 cm / 12"	Endotracheal Tube Adapters 2.5 mm, 3 mm, 3.5 mm	20
201	5 Fr	30.5 cm / 12"	Endotracheal Tube Adapters 2 mm, 2.5 mm	20
205	8 Fr	30.5 cm / 12"	Endotracheal Tube Adapters 3 mm, 3.5 mm, 4 mm	20



BALLARD™ NEONATAL AND PEDIATRIC CLOSED SUCTION SYSTEMS

MR CONDITIONAL

CLOSED SUCTION SYSTEMS - CAPPED CATHETERS

ITEM	SIZE	LENGTH	PK
20083	8 Fr	20.3 cm / 8"	20
2103	10 Fr	20.3 cm / 8"	20



MULTI-ACCESS "MAC" CATHETER

ITEM	SIZE	LENGTH	ADAPTERS	PK
1900	5 Fr	30.5 cm / 12"	Endotracheal Tube Adapters 2.5 mm, 3 mm, 3.5 mm	5
1920	5 Fr	30.5 cm / 12"	Endotracheal Tube Adapters 2.0 mm	5
1940	5 Fr	30.5 cm / 12"	Endotracheal Tube Adapters 4.0 mm	5



INTUBATION AND TRANSPORT KITS

ITEM	DESCRIPTION	PK
19000	Ballard Transport Kit for Neonates/Pediatrics, Y-Adapter, 7.5 mm	50
1902A5	Ballard Intubation Kit for Neonates/Pediatrics, Y-Adapter, 2.5 mm	50
19030	Ballard Intubation Kit for Neonates/Pediatrics, Y-Adapter, 3.0 mm	50
19035	Ballard Intubation Kit for Neonates/Pediatrics, Y-Adapter, 3.5 mm	50
1904A0	Ballard Intubation Kit for Neonates/Pediatrics, Y-Adapter, 4.0 mm	50
1904A5	Ballard Intubation Kit for Neonates/Pediatrics, Y-Adapter, 4.5 mm	50
19050	Ballard Intubation Kit for Neonates/Pediatrics, Y-Adapter, 5.0 mm	50
1906A0	Ballard Intubation Kit for Neonates/Pediatrics, Y-Adapter, 6.0 mm	50



ACCESSORIES

ITEM	DESCRIPTION	PK
11000	Ballard Manual Percussor for Neonates	12
2-020	Ballard Endotracheal Tube Adapter with Port, 2.0 mm	20
2-025	Ballard Endotracheal Tube Adapter with Port, 2.5 mm	20
2-030	Ballard Endotracheal Tube Adapter with Port, 3.0 mm	20
2-035	Ballard Endotracheal Tube Adapter with Port, 3.5 mm	20
2-040	Ballard Endotracheal Tube Adapter with Port, 4.0 mm	20
2-045	Ballard Endotracheal Tube Adapter with Port, 4.5 mm	20
1157	500 Pediatric Tracheostomy Heat and Moisture Exchanger	50

