

CAPNOFLEX™ ETCO₂ SAMPLING CANNULAS

LATEX FREE

EtCO₂ sampling cannulas for CapnoFlex™ have a proprietary design for use exclusively with the CapnoFlex module utilized with the GE Dash 3000/4000/5000™ and Carescape™ One patient monitors. They are offered in a variety of styles for all patients, in both nasal and oral-nasal designs with or without oxygen delivery. Many of the CapnoFlex kits have Nafion® sampling lines offering extended usage for patients that are receiving high humidity.



2013066-001



2013067-003



2013068-003



2013067-001VY



2013068-004

CAPNOFLEX™ LOFLO SUPPLIES, 7' TUBING

ITEM	COMPONENTS	SIZE	PK
2013066-003VY	CO ₂ Nasal Cannula	Infant	10
2013066-007	CO ₂ /O ₂ Nasal Cannula	Infant	10
2067767-006	CO ₂ /O ₂ Nasal Cannula with Nafion	Infant	10
2013066-002VY	CO ₂ Nasal Cannula	Pediatric	10
2067767-002	CO ₂ Nasal Cannula with Nafion	Pediatric	10
2013066-005	CO ₂ /O ₂ Nasal Cannula	Pediatric	10
2067767-005	CO ₂ /O ₂ Nasal Cannula with Nafion	Pediatric	10
2013067-002VY	CO ₂ Oral-Nasal Cannula	Pediatric	10
2013067-004	CO ₂ /O ₂ Oral-Nasal Cannula	Pediatric	10
2013066-001VY	CO ₂ Nasal Cannula	Adult	10
2067767-001	CO ₂ Nasal Cannula with Nafion	Adult	10
2013066-004	CO ₂ /O ₂ Nasal Cannula	Adult	10
2067767-004	CO ₂ /O ₂ Nasal Cannula with Nafion	Adult	10
2013067-001VY	CO ₂ Oral-Nasal Cannula	Adult	10
2013067-003	CO ₂ /O ₂ Oral-Nasal Cannula	Adult	10
2067767-009	CO ₂ /O ₂ Oral-Nasal Cannula with Nafion	Adult	10
2013068-001	Airway Adapter, ET > 4.0 mm	Adult/Pediatric	10
2013068-002	Airway Adapter, ET < 4.0 mm	Infant/Neonate	10
2013068-003	Airway Adapter with Nafion, ET > 4.0 mm	Adult/Pediatric	10
2013068-004	Airway Adapter with Nafion, ET < 4.0 mm	Infant/Neonate	10
2013069-001VY	Sample Line Male Luer	—	10
2013069-002VY	Sample Line with Nafion Male Luer	—	10

CAPNOFLEX™ LOFLO SUPPLIES, 12' TUBING

ITEM	COMPONENTS	SIZE	PK
1103413	CO ₂ /O ₂ Cannula*	Infant	25
1105135	CO ₂ /O ₂ Oral/Nasal Cannula*	Pediatric	25
1103408	CO ₂ Cannula*	Adult	25
1103411	CO ₂ /O ₂ Nasal Cannula*	Adult	25
1105100	CO ₂ /O ₂ Oral-Nasal Cannula*	Adult	25
1103414	Airway Adapter Set ET > 4.0 mm*	—	25
1103415	Airway Adapter Set ET < 4.0 mm*	—	25
3476-00	LoFlo Extension, 5'	—	25
1103417	Dehumidification Nafion Extension	—	10

* Must be ordered with product 3476-00—Capnoflex LoFlo extension, 5'.

CAPNOSTAT® CO₂ SENSOR FOR MAINSTREAM MONITORING

ITEM	DESCRIPTION	SIZE	PK
415036-001	Airway Adapter, Disposable, ET > 4.0 mm	Adult/Pediatric	20
415036-002	Airway Adapter, Disposable, ET < 4.0 mm	Infant/Neonate	10
412341-001	Airway Adapter, Reusable, ET > 4.0 mm	Adult/Pediatric	1

LIT-20110-C1 Rev1