

HEINE visionPRO® VIDEO LARYNGOSCOPE FOR EMS

LATEX FREE



The HEINE visionPRO® video laryngoscope is designed to provide fast, reliable, and high-quality intubation, even in challenging situations. It features an exceptionally bright and sharp image, ensuring the best view of the glottis. With its robust German engineering, the visionPRO is made to last, offering intuitive handling and easy reprocessing to ensure patient safety. Additionally, it is environmentally friendly, featuring a rechargeable battery and blades made from upcycled materials.

The visionPRO handle, display, and charger come with a 2-year warranty* for manufacturing defects.

Screen clearly displays time left and battery percentage

Screen rotates 225 degrees to reduce glare

allBRIGHT display provides true color and lifelike imaging, even in bright light

Transflective display adapts to even the brightest light source

Extra-large field of view

Removable disposable sheath in 4 configurations



MAC 3 Blade
F-000.22.253



HyMac 3 Blade
F-000.22.256

Hyperangulated blade for difficult airways

FEATURES

- The allBRIGHT display provides high-quality images with a bright and sharp video display.
- Durable construction with anodized aluminum and stainless steel components.
- Intuitive handling allows focus on intubation with quick and easy orientation.
- Components can be easily separated for cleaning and are Joint Commission compliant.
- Powered with a long-lasting, eco-friendly, and rechargeable lithium battery.
- Upcycled blades made from recycled refrigerators, significantly reducing CO₂ emissions.
- HEINE SmoothSURFACE makes the reusable handle easy-to-clean with no crevices.
- Components can be hygienically reprocessed, including STERRAD(R) sterilization for the handle.

HEINE visionPRO®

| ITEM | DESCRIPTION | PK |
|--------------|---|----|
| F-270.95.862 | visionPRO HR1 Mac Handle and visionPRO DS1 Display Unit, visionPRO CS1 Charging Station with Plug-in Power Supply | 1 |
| F-000.22.246 | Single-use HyMAC 3 Blades | 10 |
| F-000.22.256 | Single-use HyMAC 3 Blades | 50 |
| F-000.22.242 | Single-use MAC 2 Blades | 10 |
| F-000.22.252 | Single-use MAC 2 Blades | 50 |
| F-000.22.243 | Single-use MAC 3 Blades | 10 |
| F-000.22.253 | Single-use MAC 3 Blades | 50 |
| F-000.22.244 | Single-use MAC 4 Blades | 10 |
| F-000.22.254 | Single-use MAC 4 Blades | 50 |
| F-159.00.000 | Carry Bag for the visionPRO Video Laryngoscope and Accessories (Without Content) | 1 |
| F-000-22-020 | Screen Protector for Display | 5 |

HEINE visionPRO® REPLACEMENT PARTS

| ITEM | DESCRIPTION | PK |
|--------------|--|----|
| X-002.99.214 | visionPRO CS1 Charging Station with Plug-in Power Supply | 1 |
| F-007.22.350 | visionPro Display Screen DS1 | 1 |
| F-008.22.360 | visionPRO Handle HR1 | 1 |

*2-year warranty for US customers