

HEINE visionPRO® VIDEO LARYNGOSCOPE

LATEX FREE



The HEINE visionPRO® video laryngoscope is designed to provide fast, reliable, and high-quality intubation, even in challenging situations. It features an exceptionally bright and sharp image, ensuring the best view of the glottis. With its robust German engineering, the visionPRO is made to last, offering intuitive handling and easy reprocessing to ensure patient safety. Additionally, it is environmentally friendly, featuring a rechargeable battery and blades made from upcycled materials.

The visionPRO handle, display, and charger come with a 2-year warranty* for manufacturing defects.



MAC 3 Blade
F-000.22.253



HyMac 3 Blade
F-000.22.256

Hyperangulated blade for difficult airways

FEATURES

- The allBRIGHT display provides high-quality images with a bright and sharp video display.
- Durable construction with anodized aluminum and stainless steel components.
- Intuitive handling allows focus on intubation with quick and easy orientation.
- Components can be easily separated for cleaning and are Joint Commission compliant.
- Powered with a long-lasting, eco-friendly, and rechargeable lithium battery.
- Upcycled blades made from recycled refrigerators, significantly reducing CO₂ emissions.
- HEINE smoothSURFACE makes the reusable handle easy-to-clean with no crevices.
- Components can be hygienically reprocessed, including STERRAD™ sterilization for the handle.

HEINE visionPRO®

ITEM	DESCRIPTION	PK
F-270.95.862	visionPRO HR1 Mac Handle and visionPRO DS1 Display Unit, visionPRO CS1 Charging Station with Plug-in Power Supply	1
F-000.22.256	Single-use HyMAC 3 Blades	50
F-000.22.252	Single-use MAC 2 Blades	50
F-000.22.253	Single-use MAC 3 Blades	50
F-000.22.254	Single-use MAC 4 Blades	50
F-159.00.000	Carry Bag for the visionPRO Video Laryngoscope and Accessories (Without Content)	1
F-000-22-020	Screen Protector for Display	5

HEINE visionPRO® REPLACEMENT PARTS

ITEM	DESCRIPTION	PK
X-002.99.214	visionPRO CS1 Charging Station with Plug-in Power Supply	1
F-007.22.350	visionPro Display Screen DS1	1
F-008.22.360	visionPRO Handle HR1	1

*2-year warranty for US customers

LIT-20202-C1 Rev.2