

**REF** 5703

**UNITS - 1**

**RX Only**



# DISPOSABLE RESPIRATORY CIRCUITS

Adult, 72" Corrugated, Bifurcated Wye w/  
Inspiratory Ports

**CAUTION:** Use of sterilization or cleaning solutions may leave this device nonfunctional.

## SUGGESTED USE:

1. Attach the circuit pursuant to standard operating protocol.
2. Prior to patient connection:
  - a) Check that all connections are secure.
  - b) Fully extend the tubing (if using expandable circuit).
  - c) Pressurize to confirm there is no leakage or occlusion of the circuit.

**Westmed** 

5580 S. Nogales Highway  
Tucson, Arizona 85706 USA  
Phone: 800-975-7987  
Fax: 520-294-6061  
www.westmedinc.com

Assembled in Mexico  
LABEL PN 77090, Rev. 04



**LOT**



MT Promedt Consulting GmbH  
Altenhofstr. 80  
66386 St. Ingbert, Germany  
Patent Pending