

REF 7-7105-66


UNITS - 1
RX Only



C **€**0482

DISPOSABLE ANESTHESIA BREATHING CIRCUITS

Jackson Rees Circuit, 12", Norman Mask
Elbow, 1/2L Not Made with Latex Bag,
Bleed Tail Valve, Fresh Gas Line

 **CAUTION:** Use of sterilization or cleaning solutions may leave this device nonfunctional.

SUGGESTED USE:

1. Attach the circuit pursuant to standard operating protocol.
2. Prior to patient connection: Check that all connections are secure. Fully extend the tubing (if using expandable circuit). Pressurize to confirm there is no leakage or occlusion of the circuit.

Westmed 

5580 S. Nogales Highway
Tucson, Arizona 85706 USA
Phone: 800-975-7987
Fax: 520-294-6061
www.westmedinc.com



MT Promedt Consulting GmbH
Altenhofstr. 80
66386 St. Ingbert, Germany
Patent Pending

Assembled in Mexico
LABEL PN 79000, Rev. 01



LOT