



5
Adult

Top Valve Vent Mask™

Premium Inflatable Non-Conductive Anesthesia Mask

REF VR5500 / 5 / Adult

Directions for Use

Inspect cushion for proper inflation pressure before use. Use luer tip syringe to adjust cushion pressure, if needed.



Caution

Do not attempt to sterilize, autoclave or clean this device.



R_X ONLY



Ventlab, LLC
2710 Northridge Dr. NW, Suite A
Grand Rapids, MI 49544 USA



Manufactured for Ventlab, LLC,
a Subsidiary of SunMed Holdings, LLC
www.Sun-Med.com
Made in China by SunMed AP, LTD

VR5500.IFU Rev.3 2022/11