

## LIMB-O® SINGLE LIMB WARM BREATHING CIRCUIT

LATEX FREE | DEHP FREE

Limb-O® is a single divided tube utilizing a flexible septum running the length of the tube. The septum effectively transfers heat from the expiratory to the inspiratory side of the circuit. Increased thermal efficiency rates during surgery keeps gases and supports normothermic efforts.



### FUNCTIONALITY

- The septum in Limb-O decreases compliance in the circuit to significantly increase delivered volume without sacrificing drape
- Effectively transfers heat from the expiratory to the inspiratory side of the circuit
- Increased thermal efficiency rates during surgery keeps gases warm
- Less weight means less torque on the ET tube and reduced plastic waste
- Accommodates both adult and pediatric patients



Single limb circuit divided by a flexible septum

# LIMB-O™ SINGLE LIMB WARM BREATHING CIRCUIT

LATEX FREE | DEHP FREE

## LIMB-O™ SINGLE LIMB BREATHING CIRCUIT

ITEM	SIZE	COMPONENTS	PK
AFN119X4	72" Circuit	3L, Sampling Elbow, HEPA, Size 5 Mask, & M/M GSL	20
AFN1292X	72" Circuit	3L, Sampling Elbow, HEPA, Size 5 Mask, & HME	20
AFN510X4	72" Circuit	3L, Sampling Elbow, 2 B/V Filters, Size 5 Mask, & M/M GSL	20
AFN52004	72" Circuit	3L, Sampling Elbow, 2 B/V Filters, Size 5 Mask, HMEF, & M/M GSL	20
AFN529X4	72" Circuit	3L, Sampling Elbow, 2 B/V Filters, Size 5 Mask, & M/M GSL	20
AFN560X4	72" Circuit	3L, Sampling Elbow, 2 B/V Filters, Size 5 Mask, & M/M GSL	20
AFN580AC	72" Circuit	2L, Sampling Elbow, 2 B/V Filters, Size 5 Mask, HME, & Integrated GSL	20
AFN5QXX4	72" Circuit	1L, Sampling Elbow, 2 B/V Filters, & M/M GSL	20
AFN5X0TZXX6	72" Circuit	3L, Sampling Elbow, 2 B/V Filters, Size 5 Mask, HME, & M/M GSL	20
AFND8XX4P	72" Circuit	2L, Sampling Elbow, B/V, HMEF, & M/M GSL	20
AFNX1900	72" Circuit	3L, Sampling Elbow, Size 5 Mask, HMEF, & M/M GSL	20
AFNX2004	72" Circuit	3L, Sampling Elbow, Size 5 Mask, HMEF, & M/M GSL	20
AFNX2044	72" Circuit	3L, Sampling Elbow, Size 5 Mask, HMEF, & M/M GSL	20
AFNX20X4	72" Circuit	3L, Sampling Elbow, Size 5 Mask, & M/M GSL	20
AFNX20XZP	72" Circuit	3L, Sampling Elbow, Size 5 Mask, HMEF, & M/M GSL	20
AFNX2904	72" Circuit	3L, Sampling Elbow, Size 5 Mask, HMEF, & M/M GSL	20
AFNX8000	72" Circuit	2L, Sampling Elbow, Size 5 Mask, HMEF, & M/M GSL	20
AFNX8044	72" Circuit	2L, Sampling Elbow, Size 5 Mask, HMEF, & M/M GSL	20
AFNXOXXX	72" Circuit	2L, Sampling Elbow	20
AFNXX004XX6	72" Circuit	3L, Sampling Elbow, Size 5 Mask, HMEF, M/M GSL, & ABS bag arm	20
AFNXXXXX	72" Circuit	Sampling Elbow	20
AFPXXXXX	108" Circuit	Sampling Elbow	20
AGNX204Z	72" Circuit	3L, Plain Elbow, Size 5 Mask, HMEF, & M/M GSL	15
AGNX2944	72" Circuit	3L, Plain Elbow, Size 5 Mask, HMEF, & M/M GSL	20
AGNXOXXC	72" Circuit	2L, Plain Elbow, & Integrated GSL	20
AGPX2044	108" Circuit	3L, Plain Elbow, Size 5 Mask, HMEF, & M/M GSL	20
AGPX2XXX	108" Circuit	3L, Plain Elbow	20
ALY320X4	72" Circuit	3L, Sampling Elbow, HEPA, Size 5 Mask, & M/M GSL	20
AMN520X4	72" Circuit	3L, Sampling Elbow, 2 B/V Filters, Size 5 Mask, & M/M GSL	20
AMN520X4F	72" Circuit	3L, Sampling Elbow, 2 B/V Filters, Size 5 Mask, HMEF, & M/M GSL	20
AMN580X4	72" Circuit	2L, Sampling Elbow, 2 B/V Filters, Size 5 Mask, & M/M GSL	20
AXXXXX61	72" Circuit	1L, Plain Elbow, & Integrated GSL	20
G15LXXXX	72" Circuit	3L, Sampling Elbow, Size 5 Mask, HMEF, M/M GSL, & Suction Tubing	20
GASFXXXX	72" Circuit	3L, Sampling Elbow, 2 B/V Filters, HME, M/M GSL, Yankauer, & Suction Tubing	20
GD6FXXXX	72" Circuit	3L, Sampling Elbow, Size 5 Mask, HMEF, M/M GSL, Yankauer, & Suction Tubing	20