

RADIOLUCENT LEADWIRES WITH MULTI-LINK™ X2 SYSTEM AHA

SINGLE-SOURCE SOLUTION FOR ALL ECG NEEDS

Unlike many other ECG solution providers, we offer not only reusable and single-patient use portfolios, but also radiolucent leadwires. By standardizing to Multi-Link™ X2, you can choose the Multi-Link leadwires that fit your needs throughout your facility. The Multi-Link X2 ECG System is your single-source solution for all ECG needs.

RADIOLUCENT LEADWIRES FOR X-RAY

As part of the standardized solution, you can use standard reusable Multi-Link X2 ECG cables with the Multi-Link radiolucent leadwires. Choosing your ECG solution for X-ray environment is easy:



SELECT YOUR 3.6 M/12 FT CABLE:

MONITOR INPUT	SKU	MONITOR BRAND	LEADS	TYPE
	2017003-001	GE HealthCare	3 or 5	Cable
	2017005-001		6	Cable
	2017006-001		12	Cable
	2090601-013	GE HealthCare CARESCAPE ONE	12	Cable
	2090100-052	Philips	3 or 5	Cable
	2090100-054		6	Cable
	2090100-056		12	Cable
	2090101-091	Dräger	5 or 6	Adapter
	2090100-062	Mindray Traditional	3 or 5	Cable
	2090100-066		12	Cable
	2090600-012	Mindray Passport V12/21	5	Cable
	2090600-022	Mindray Datascope	5	Cable
	2090100-072	Spacelabs	5	Cable
	2090100-076		12	Cable
	2090600-038	Physio Control	5	Cable
	2090600-062	Zoll R-Series	5	Cable
	2090600-052	Zoll X-Series	5	Cable

KEY FEATURES

- Multi-Link X2 compatible
- Free of phthalates such as DEHP*
- Does not contain BPA*
- Latex Free
- Designed by ECG experts
- Tested to meet/exceed ECG standards
- Compatibility tested with systems

* These products are not manufactured with natural rubber latex and do not contain DEHP, BBP, DBP or DIBP in concentrations above 0.1% w/w

ADD YOUR RADIOLUCENT LEADWIRE:

	SKU	LEADS	LENGTH
	403751-023	5	150 cm/59 in
	403751-024	10	