

TRI-FLO® SUCTION CATHETERS

LATEX FREE

Tri-Flo catheters are designed with two triangular eyes and a beveled tip to increase suctioning efficiency while reducing patient trauma.

Tri-Flo catheters are available in looped and straight packaging configurations.




T60C




Beveled tip design to increase suctioning efficiency



Top - Control Port
Bottom - Without Control Port
For continuous or controlled suctioning



T64C



Depth markings from 5 to 20 cm in 2 cm increments are standard on all 5-6, 8 and 10 Fr sizes to help facilitate proper catheter placement.

Clear connectors are standard on all 5/6 and 8 Fr sizes to help easily visualize patient secretions.

TRI-FLO® SINGLE CATHETERS

ITEM	DESCRIPTION	SIZE	PK
T61	Looped Catheter Without Control Port	10 Fr	50
T60	Looped Catheter Without Control Port	14 Fr	50
T62	Looped Catheter Without Control Port	18 Fr	50
T63C	Looped Catheter With Control Port	5-6 Fr	50
T64C	Looped Catheter With Control Port	8 Fr	50
T61C	Looped Catheter With Control Port	10 Fr	50
T68C	Looped Catheter With Control Port	12 Fr	50
T60C	Looped Catheter With Control Port	14 Fr	50
T62C	Looped Catheter With Control Port	18 Fr	50
T260	Straight Catheter Without Control Port	14 Fr	100
T263C	Straight Catheter With Control Port	5-6 Fr	100
T264C	Straight Catheter With Control Port	8 Fr	100
T261C	Straight Catheter With Control Port	10 Fr	100
T268C	Straight Catheter With Control Port	12 Fr	100
T260C	Straight Catheter With Control Port	14 Fr	100
T262C	Straight Catheter With Control Port	18 Fr	100
T164C	No-Touch Catheter With Pop-Up Basin	8 Fr	100
T161C	No-Touch Catheter With Pop-Up Basin	10 Fr	100
T168C	No-Touch Catheter With Pop-Up Basin	12 Fr	100
T160C	No-Touch Catheter With Pop-Up Basin	14 Fr	100
T162C	No-Touch Catheter With Pop-Up Basin	18 Fr	100

NO-TOUCH CATHETER



T160C

Packaged in a plastic sleeve to add a layer of safety between the clinician and the patient.

Pop-up basin for your solution.

TRI-FLO® CATH-N-GLOVE® KITS

Economy Cath-N-Glove® suction kits include a Tri-Flo® suction catheter, a pop-up solution basin, and one or two powder-free vinyl exam gloves. These kits are conveniently enclosed in a plastic peel pouch for quick and easy preparation. Wallet Catheter kits include a Tri-Flo suction catheter, a pop-up solution basin, and two precuffed exam gloves, neatly arranged in an inner wallet.



4664A

TRI-FLO® CATH-N-GLOVE™ CATHETER KITS

ITEM	DESCRIPTION	SIZE	PK
4863T	Wallet Catheter Kit Two Nitrile Gloves	5-6 Fr	100
4867T	Wallet Catheter Kit Two Nitrile Gloves	8 Fr	100
4865T	Wallet Catheter Kit Two Nitrile Gloves	10 Fr	100
4864T	Wallet Catheter Kit Two Nitrile Gloves	14 Fr	100
4868T	Wallet Catheter Kit Two Nitrile Gloves	12 Fr	100
4893T	Economy Catheter Kit Two Vinyl Gloves	5-6 Fr	100
4897T	Economy Catheter Kit Two Vinyl Gloves	8 Fr	100
4895T	Economy Catheter Kit Two Vinyl Gloves	10 Fr	100
4898T	Economy Catheter Kit Two Vinyl Gloves	12 Fr	100
4894T	Economy Catheter Kit Two Vinyl Gloves	14 Fr	100
4896T	Economy Catheter Kit Two Vinyl Gloves	18 Fr	100
4697T	Economy Catheter Kit 1 Vinyl Glove	8 Fr	100
4695T	Economy Catheter Kit 1 Vinyl Glove	10 Fr	100
4698T	Economy Catheter Kit 1 Vinyl Glove	12 Fr	100
4694T	Economy Catheter Kit 1 Vinyl Glove	14 Fr	100