

# ADJUSTABLE PEEP VALVE

LATEX FREE | SINGLE USE | MR CONDITIONAL

- Designed for use with AirLife™ AirFlow manual resuscitators, manual resuscitators, and ventilator circuits with a 30 mm exhalation port.
- Pressure range: 5-20 cm H<sub>2</sub>O
- Large print on scale for improved visibility of setting
- Connects to exhalation port of manual resuscitator
- MR Conditional up to 7 tesla



VP700



VP702

ITEM	DESCRIPTION	SIZE	PK
VP700	Adjustable PEEP Valve, 5-20 cm H <sub>2</sub> O	19 mm	20
VP702	Adjustable PEEP Valve + Adapter, Converts to 30 mm	19 mm + 30 mm Adapter	20