



RESPIRGARD™ NEBULIZER SYSTEM

LATEX FREE | SINGLE PATIENT USE | DISPOSABLE

The AirLife™ filtered nebulizer system was designed and tested for the clinical investigation of aerosolized pentamidine; however, it is suitable to deliver most any kind of aerosolized respiratory medication. Its design and configuration help minimize the environmental impact of any aerosolized medication compared to a standard small volume nebulizer.



ITEM	DESCRIPTION	PK
124030EU	AirLife filtered nebulizer system	20
PERFORMANCE SPECS		
MMAD	3.5 +/- 0.2 μ	
% Respirable Mass (1-5 μ)	54 +/- 0.4%	
Respirable Mass Output rate (1-5 μ)	83 +/- 8 μ g/min	
Treatment Time	7.3 +/- 0.4 mins @ sputter + 1 min	
Residual Volume	1.0 +/- 0.1 g	