


SUNONE™ SYSTEM GREENLINE®/D™ FIBER OPTIC GREEN HEEL LARYNGOSCOPE BLADES

LATEX FREE | STERILE | SINGLE PATIENT USE | DISPOSABLE

GREENLINE®/D™ GREEN HEEL BLADES





GREEN HEEL MACINTOSH

ITEM	DESCRIPTION	LENGTH	HEIGHT	SIZE
5-5332-01	Infant	92 mm	22 mm	1
5-5332-02	Child	100 mm	21 mm	2
5-5332-03	Medium Adult	130 mm	22 mm	3
5-5332-35	Ext. Med. Adult	144 mm	22 mm	3.5
5-5332-04	Large Adult	155 mm	25 mm	4

GREEN HEEL MILLER

ITEM	DESCRIPTION	LENGTH	HEIGHT	SIZE
5-5333-20	Premature	67 mm	10 mm	00
5-5333-00	Neonate	77 mm	10 mm	0
5-5333-01	Infant	102 mm	10 mm	1
5-5333-02	Child	153 mm	12 mm	2
5-5333-03	Medium Adult	195 mm	14 mm	3
5-5333-04	Large Adult	205 mm	18 mm	4



FLEX TIP BLADES

ITEM	DESCRIPTION	LENGTH	HEIGHT	SIZE
5-5382-03	Medium Adult	195 mm	18 mm	3
5-5382-04	Large Adult	205 mm	18 mm	4

ROBERTSHAW BLADES

ITEM	DESCRIPTION	LENGTH	HEIGHT	SIZE
5-5337-00	Neonate	94 mm	11 mm	0
5-5337-01	Infant	112 mm	13 mm	1