


TRACH-VAC SUBGLOTTIC SUCTION ENDOTRACHEAL TUBE WITH FLEX-TIP™

LATEX FREE | SINGLE PATIENT USE | DISPOSABLE

Incorporate the ease of Flex-Tip™ to your most challenging ICU intubation events and reduce intubation time and insertional trauma for your most vulnerable patients. 6-28% of ICU intubations are rated as "difficult." The subglottic suction of the Flex-Tip Endotracheal Tube is proven to decrease VAP rates in ventilated patients¹, making for less traumatic insertions and increases "first pass success"² rates.



FLEX-TIP™
Endotracheal Tube Tip

STANDARD
Endotracheal Tube Tip

ITEM	DESCRIPTION	CUFF DIA.	ID	OD	PK
I-SFTVVC-60-10	TrachVac with PVC Cuff	25 mm	6.0 mm	9.0 mm	10
I-SFTVVC-65-10	TrachVac with PVC Cuff	25 mm	6.5 mm	9.8 mm	10
I-SFTVVC-70-10	TrachVac with PVC Cuff	26 mm	7.0 mm	10.4 mm	10
I-SFTVVC-75-10	TrachVac with PVC Cuff	26 mm	7.5 mm	11.2 mm	10
I-SFTVVC-80-10	TrachVac with PVC Cuff	28 mm	8.0 mm	11.8 mm	10
I-SFTVVC-85-10	TrachVac with PVC Cuff	28 mm	8.5 mm	12.6 mm	10
I-SFTVVC-90-10	TrachVac with PVC Cuff	28 mm	9.0 mm	13.1 mm	10
I-SFTVPU-60-10	TrachVac with Polyurethane Cuff	20 mm	6.0 mm	9.0 mm	10
I-SFTVPU-65-10	TrachVac with Polyurethane Cuff	24 mm	6.5 mm	9.8 mm	10
I-SFTVPU-70-10	TrachVac with Polyurethane Cuff	24 mm	7.0 mm	10.4 mm	10
I-SFTVPU-75-10	TrachVac with Polyurethane Cuff	26 mm	7.5 mm	11.2 mm	10
I-SFTVPU-80-10	TrachVac with Polyurethane Cuff	26 mm	8.0 mm	11.8 mm	10
I-SFTVPU-85-10	TrachVac with Polyurethane Cuff	30 mm	8.5 mm	12.6 mm	10
I-SFTVPU-90-10	TrachVac with Polyurethane Cuff	30 mm	9.0 mm	13.1 mm	10