

# HYPERINFLATION SYSTEM

LATEX FREE | SINGLE PATIENT USE | DISPOSABLE

The Hyperinflation System is available in four different pleated inflation bag sizes. Available with the Ventlab™ signature integrated pressure manometer and components to create a custom system.



## HYPERINFLATION SYSTEM

- System equipped with a stay-put dial to enable clinician to achieve targeted static pressure
- Four latex-free pleated breathing bag sizes available: 1/2, 1, 2, 3 liter
- Seven inflatable face mask sizes available: small infant, infant, small child, child, adult (regular, large, wide)
- Non-inflatable face masks can be specified

ITEM	DESCRIPTION	PK
HS1000 Series	Hyperinflation System	20
HS4000 Series	Hyperinflation System + Integrated Manometer	20

LIT-20021-C1 Rev.1