

VITAL SIGNS® PRESSURE INFUSOR

LATEX FREE | SINGLE PATIENT USE

The durable construction of the Vital Signs® pressure infusors helps maintain constant pressure and prevents leaks. It also comes with a nylon netting on the back ensures simple and secure placement of the fluid bag.



FUNCTIONALITY

- Robust and durable construction to hold constant pressure and prevent leaks
- Single patient use to minimize the spread of nosocomial infections via touch contamination
- Color-coded printing for rapid identification of size
- Materials do not contain natural rubber latex or neoprene
- Available in 3 sizes – 500 mL, 1000 mL and 3000 mL

PRESSURE INFUSOR

ITEM	DESCRIPTION	SIZE	PK
IN800012	Pressure Infusor, Netted Backing	500 mL	12
IN900012	Pressure Infusor, Netted Backing	1000 mL	12
IN950012	Pressure Infusor, Netted Backing	3000 mL	12