

ANESTHESIA CIRCUITS



Dual Limb Anesthesia Circuits

LATEX FREE

- Configurations available with face masks, filters, heat moisture exchangers with filters (HMEFs), and gas sampling lines
- Designed and tested to ensure accurate readings and peak performance
- ISO-certified, available in neonatal, pediatric, and adult sizes
- Expandable dual-limb hoses expand to multiple lengths
- 100% leak tested



ITEM	DESCRIPTION	SIZE	BAG	LENGTH	PK
A4UX2044	Expandable Anesthesia Circuit, Parallel Wye, Sampling Elbow, Size 5 Mask, HMEF, and M/M GSL	Adult	3 Liter	108"	20
ADU520X4D	Expandable Anesthesia Circuit, Ported Parallel Wye, Plain Elbow, 2 B/V Filters, Size 5 Mask, HMEF, and M/M GSL	Adult	3 Liter	108"	20
AMY520X4F	Expandable Anesthesia Circuit, Parallel Wye, Sampling Elbow, 2 B/V Filters, Size 5 Mask, HMEF, and M/M GSL	Adult	3 Liter	120"	20
A4UX6044	Expandable Anesthesia Circuit, Parallel Wye, Sampling Elbow, Size 5 Mask, HMEF, and M/M GSL	Adult	3 Liter	108"	20
A5U510X4	Expandable Anesthesia Circuit, Ported Parallel Wye, Sampling Elbow, 2 B/V Filters, Size 5 Mask, and M/M GSL	Adult	3 Liter	108"	20
A4FX2004	Expandable Anesthesia Circuit, Parallel Wye, Sampling Elbow, Size 5 Mask, HMEF, and M/M GSL	Adult	3 Liter	75"	20
A5Z11XXX	Expandable Anesthesia Circuit, Ported Parallel Wye, Sampling Elbow, and B/V	Adult	3 Liter	90"	20
A5Z329X4	Expandable Anesthesia Circuit, Ported Parallel Wye, Sampling Elbow, B/V, Size 5 Mask, and M/M GSL	Adult	3 Liter	90"	20
A5Z12XXX	Expandable Anesthesia Circuit, Ported Parallel Wye, Sampling Elbow, and B/V	Adult	3 Liter	96"	20
B17B4XX4	Expandable Anesthesia Circuit, Parallel Pediatric Wye, Sampling Elbow, B/V, and M/M GSL	Pediatric	1 Liter	108"	20
1552017	Expandable Anesthesia Circuit, Ported Neonatal Wye, Sampling Elbow, B/V, HMEF, and M/M GSL	Neonate	1/2 Liter	60"	20
B1314XXX	Corrugated Anesthesia Circuit, Parallel Pediatric Wye, Sampling Elbow, and B/V	Pediatric	1 Liter	60"	20



Single Limb Anesthesia Circuits

Uni-Lim™ and Pedi-Lim™ Single Limb Anesthesia Circuits

LATEX FREE | SINGLE PATIENT MULTIPLE USE

The Uni-Lim™ and Pedi-Lim™ Single Limb Coaxial Anesthesia Circuit features a tube-within-a-tube design that works to eliminate twisting between the inspiratory and expiratory limb. It provides natural heat and moisture conservation properties.

- Configurations available with face masks, filters, heat moisture exchangers with filters (HMEFs), and gas sampling lines



9810 Uni-Lim



9812 Pedi-Lim



M/M 10' Gas Sampling Line

Breathing bag options in 1, 2, and 3 L



M/M 10' Gas Sampling Line

Breathing bag options in 1, 2, and 3 L

ITEM	DESCRIPTION	SIZE	BAG	LENGTH	FILTER	PK
9812	Pedi-Lim Anesthesia Circuit, Sampling Elbow, and M/M GSL	Pediatric	1 Liter	60"	2 B/V	20
5610004	Uni-Lim Anesthesia Circuit, Sampling Elbow, Size 5 Mask, HMEF with Port, and M/M GSL	Adult	3 Liter	72"	2 B/V	20
9694	Uni-Lim Anesthesia Circuit, Sampling Elbow, and M/M GSL	Adult	3 Liter	72"	2 B/V	20
9810	Uni-Lim Anesthesia Circuit, Sampling Elbow, and M/M GSL	Adult	3 Liter	72"	2 B/V	20

LIMB-O® Single Limb Circuits

LATEX FREE

The LIMB-O® Single Limb Circuit is divided by a flexible septum. It effectively transfers heat from the expiratory to the inspiratory side of the circuit.

- Increased thermal efficiency rates during surgery keeps gases warm and body temperature normothermic
- Less weight means less torque on the ETT and less plastic waste
- Decreased compliance helps lead to accurate delivered volume
- LIMB-O circuit services patients from pediatric to adult
- Configurations available with face masks, filters, heat moisture exchangers with filters (HMEFs), and gas sampling lines



ITEM	DESCRIPTION	BAG	LENGTH	PK
AFNX2044	Anesthesia Circuit, LIMB-O, Sampling Elbow, Size 5 Mask, HMEF, & M/M GSL	3 Liter	72"	20
AMN520X4F	Anesthesia Circuit, LIMB-O, Sampling Elbow, 2 B/V Filters, Size 5 Mask, HMEF, & M/M GSL	3 Liter	72"	20

Jackson-Rees Anesthesia Circuits

The Jackson-Rees circuit, also referred to as a Mapleson F circuit, provides a direct connection between the fresh gas hose and fresh gas outlet on the anesthesia machine.



ITEM	DESCRIPTION	SIZE	BAG	LENGTH	PK
5101ELW1	Anesthesia Circuit, 10" Jackson-Rees, 1L, Sampling Elbow, Washington Tee, & O ₂ Tubing	Adult	1 Liter	11"	20

Circuit Accessories

Breathing Bags

LATEX FREE | SINGLE PATIENT MULTIPLE USE

ITEM	DESCRIPTION	SIZE	PK
5060NL	Breathing Bag	1/2 Liter	40
5061NL	Breathing Bag	1 Liter	40
5062NL	Breathing Bag	2 Liter	40
5063NL	Breathing Bag	3 Liter	40



5060NL

5061NL

5062NL

5063NL

Anesthesia Accessories



FH106347



FH106348



FH106380



Circuit Extension Sets

DISPOSABLE | INDIVIDUAL PACKAGING
LATEX FREE | SINGLE PATIENT MULTIPLE USE

ITEM	DESCRIPTION	LENGTH	PK
5108EX	Anesthesia Circuit, Expandable Extension	108"	20
5024EX	Anesthesia Circuit, Expandable Extension	24"	20

Flex Tubing (Catheter Mounts)

DISPOSABLE | INDIVIDUAL PACKAGING
LATEX FREE | SINGLE PATIENT MULTIPLE USE

Flex tubing for HME/HMEF is available for use with hydroscopic condenser humidifiers. Made of disposable, latex free material.

- For use with hygroscopic condenser humidifiers
- Standard 6" x 15 mm ID flex tubing

ITEM	DESCRIPTION	SIZE	PK
FH106347	6" x 15 mm ID Flex Tubing	15 mm OD x 15 mm OD	50
FH106348	6" x 15 mm ID Flex Tubing	15 mm OD x 22 mm OD/15 mm ID	50
FH106380	6" x 15 mm ID Flex Tubing	15 mm OD x 22 mm OD/15 mm ID 15 mm OD Connector	50




Straight T Airway Adapters

The Straight T Airway Adapter allows for gas sampling while airway procedures are being performed.

ITEM	DESCRIPTION	PK
73385-HEL	Vital Signs Straight T Adapter with Female Luer Lock Connector, Disposable	50

Vent Mask™ Inflatable Anesthesia Masks

- Superior ergonomic cushion design for an optimal seal with minimum pressure
- Flexible body allows gentle deep seating of cushion to face
- Inflation valve adjusts cuff pressure for optimal seal
- Ample nose bridge room for proper fit
- Centered connector port for secure gripping and easy connect/disconnect
- Removable hook ring enables manual or head strap application
- Premium soft cushioned mask
- Low-pressure soft seal with superior seating characteristics



ITEM	DESCRIPTION	SIZE	VALVE	HOOK RING	PK
VR1100-20	Small Infant Vent Mask	Size 1	Tail	—	20
VR2100-20	Infant Vent Mask	Size 2	Tail	—	20
VR3100-20	Small Child Vent Mask	Size 3	Tail	Grey	20
VR4100-20	Child Vent Mask	Size 4	Tail	Grey	20
VR5000-20	Wide Regular Adult Vent Mask	Size 5W	Tail	Grey	20
VR5100-20	Regular Adult Vent Mask	Size 5	Tail	Grey	20
VR6100-20	Adult Large Vent Mask	Size 6	Tail	Grey	20

Vent Mask II™ Anesthesia Masks

- Non-conductive anesthesia mask
- Air-filled soft cushioned mask
- Pliable ultra-thin wall
- Low-pressure soft seal with superior seating characteristics



ITEM	DESCRIPTION	SIZE	VALVE	HOOK RING	PK
VR0040-20	Extra Small Neonate Blow Mold	Size 00	—	—	20
VR0400-20	Small Neonate Blow Mold Mask	Size 0	—	—	20
VR1400-20	Small Infant Blow Mold Mask	Size 1	—	—	20
VR2400-20	Infant Blow Mold Mask	Size 2	—	—	20
VR3400-20	Small Child Blow Mold Mask	Size 3	—	Grey	20
VR4400-20	Child Blow Mold Mask	Size 4	—	Grey	20
VR5400-20	Adult Blow Mold Mask	Size 5	—	Grey	20
VR6400-20	Large Adult Blow Mold Mask	Size 6	—	Grey	20
VR7400-20	Extra Large Wide Adult Blow Mold Mask	Size 5W	—	Grey	20

Flexible Anesthesia Masks

The Flexible Anesthesia Mask comes packaged in many different colors to ensure quick and easy size identification. The teardrop-shaped mask features a cushion that is easy to deflate or inflate.

- Easy to inflate or deflate air cushion with top valve convenience
- Soft face plate
- Superior grip
- Pediatric sizes offered in bubble gum scent



ITEM	DESCRIPTION	SIZE	PK
6810EU	Flex Dome, Top Valve, Air Filled Cushion, Clear	1 - Neonate	20
6820	Flex Dome, Top Valve, Air Filled Cushion, Clear	2 - Infant	20
6821	Flex Dome, Top Valve, Air Filled Cushion, Clear	2 - Infant	100
6830	Flex Dome, Top Valve, Air Filled Cushion, Clear	3 - Toddler	20
6831	Flex Dome, Top Valve, Air Filled Cushion, Clear	3 - Toddler	100
6840	Flex Dome, Top Valve, Air Filled Cushion, Clear	4 - Child	20
6841	Flex Dome, Top Valve, Air Filled Cushion, Clear	4 - Child	100
6850	Flex Dome, Top Valve, Air Filled Cushion, Clear	5 - Adult	20
6851	Flex Dome, Top Valve, Air Filled Cushion, Clear	5 - Adult	100
6860	Flex Dome, Top Valve, Air Filled Cushion, Clear	6 - XL Adult	20
6861	Flex Dome, Top Valve, Air Filled Cushion, Clear	6 - XL Adult	100

Items 6820, 6821, 6830, 6831, 6840, and 6841 are also available bubble gum scented.

CO₂ Absorbents

Medisorb® Medical Soda Lime CO₂ Absorbent

LATEX FREE | DISPOSABLE | DEHP FREE

The Medisorb® CO₂ Absorbent has a unique granular shape that allows better gas flow through the absorbent bed. This product comes in a variety of packaging options for every clinical need.

- Medisorb has been specified for use with and validated for optimal performance with GE Healthcare® anesthesia machines
- Medisorb EX and low-alkaline Medisorb EF EX disposable canisters provide easy installation and quick change mid-procedure for optimal utilization
- Medisorb EF is free from a strong alkali, meaning it does not generate Compound A or carbon monoxide
- Medisorb EF offers a permanent color change to violet when desicated



M1173310

2079796-001

ITEM	DESCRIPTION	PK
M1173310	Medisorb Original (for use with Aespire, Avance, Aisys and ADU/Stien)	6
M1173311	Medisorb EF (for use with Aespire, Avance, Aisys and ADU/Stien)	6
2079796-001	Medisorb EX (for use with Carestation 600 Series)	6
2079797-001	Medisorb EF EX (for use with Carestation 600 Series)	6
427000100	Medisorb Pre-Packed Rounds, Disposable	12
8570043	Medisorb Twin Pack	2
1407-3201-000	Dust Filter	40



BRO-20085-C1 Rev. 1

