


BALLARD™ 7DAY CLOSED SUCTION SYSTEM

STERILE | SINGLE PATIENT USE

Introducing the newest addition to the Ballard Closed Suction System portfolio. The Ballard 7Day Closed Suction System is a simplified, streamlined, easier-to-use portfolio of products for patients with an artificial airway for more than 72 hours. With a new look, clear adapter, and turbulent cleaning technology, these products empower clinicians with the flexibility to change out the catheter as needed, up to 7 days use, while retaining confidence in the quality and reliability of the product they are accustomed to with the Ballard brand. **MR safe and not made with latex/DEHP.**



Available in T-Piece Styles
Tracheostomy Length

- Tethered T-Piece cap
- Numerical indices on all sizes

Double Swivel Elbow Style
Endotracheal Length

ITEM	DESCRIPTION	DURATION	SIZE	ADAPTER TYPE	ET/TRACH	LENGTH	ACCESSORIES	PK
727101	Ballard 7Day Closed Suction System	7 Days	10 Fr	Double Swivel Elbow	ET Length	21.3" / 54 cm	15 mm x 22 mm Flex Connector	20
727121	Ballard 7Day Closed Suction System	7 Days	12 Fr	Double Swivel Elbow	ET Length	21.3" / 54 cm	15 mm x 22 mm Flex Connector	20
727141	Ballard 7Day Closed Suction System	7 Days	14 Fr	Double Swivel Elbow	ET Length	21.3" / 54 cm	15 mm x 22 mm Flex Connector	20
727161	Ballard 7Day Closed Suction System	7 Days	16 Fr	Double Swivel Elbow	ET Length	21.3" / 54 cm	15 mm x 22 mm Flex Connector	20
7271035	Ballard 7Day Closed Suction System	7 Days	10 Fr	T-Piece	Trach Length	12" / 30.5 cm	22 mm Flex Tube	20
7271235	Ballard 7Day Closed Suction System	7 Days	12 Fr	T-Piece	Trach Length	12" / 30.5 cm	22 mm Flex Tube	20
7271435	Ballard 7Day Closed Suction System	7 Days	14 Fr	T-Piece	Trach Length	12" / 30.5 cm	22 mm Flex Tube	20
7271635	Ballard 7Day Closed Suction System	7 Days	16 Fr	T-Piece	Trach Length	12" / 30.5 cm	22 mm Flex Tube	20