


HCA FORMULARY REORDER GUIDE

ECG MONITORING



HCA has entered into a formulary for bedside and telemetry ECG monitoring with AirLife® and their Multi-Link™ X2 ECG System. Below is the best direct cross for HCA facilities based on historical usage.

Current Product			AirLife Product							Cross Reference	
Manufacturer	Item	Description	Item	Description	Product Category	Reusable/ Disposable	UOM	QTY	Cross Type	Comments	
Cardinal Health 200 LLC	33526T	ADAPTER EKG 2 CONN PHLPS INTEL TELEM	2090101-050	Adapter, ECG, Philips Intellivue 5/6-Lead	Adapter, ECG	Reusable	Each	1	Clinically Equivalent	Original is adapter, we offer cables	
Cardinal Health 200 LLC	33507	ADAPTER EKG 12 PIN KDL DL LF 5 LD MNTR	2090100-052	Cable, ECG, Philips 3/5-Lead, AHA, 3.6 m/12'	Cable, ECG	Reusable	Each	1	Clinically Equivalent	Original is adapter, we offer cables	
Cardinal Health 200 LLC	33523T	ADAPTER SWITCHING KENDALL LD ECG NONSTRL REUSABLE DL 5-LEAD	—	—	—	—	—	—	No Cross	No adapters needed for GE ApexPro CH	
Cardinal Health 200 LLC	33532T	ADAPTER EKG PRO NON-INVASIVE CABLE ELECT PHYS MEDICINE KDL	—	—	—	—	—	—	No Cross	No adapters needed for GE ApexPro CH	
ConMed Corp.	38-406-5	CABLE ECG NON-INVASIVE CV NONSTRL DISP	2017003-001	Cable, ECG, GE, 3/5-Lead, AHA, 3.6 m/12'	Cable, ECG	Reusable	Each	1	Clinically Equivalent	Original is for ConMed Backpad, offering cable for standard leadwires	
Curbell Medical Products Inc.	CBT-10NR-10CM-0001	CBL 10LD GE MULTILINK STANDARD TRUNK	2017003-001	Cable, ECG, GE, 3/5-Lead, AHA, 3.6 m/12'	Cable, ECG	Reusable	Each	1	Exact	—	
Curbell Medical Products Inc.	TCB-03NA-10AS-0010	CBL 3LD GE TO SAFETY DIN STANDARD TRUNK	2017004-001	Cable, ECG, Nihon Kohden, GE, 3-Lead, AHA, 3.6 m/12'	Cable, ECG	Reusable	Each	1	Exact	—	
Curbell Medical Products Inc.	TCB-05ND-10AS-0000	CABLE EKG GRN 5 LD 11 PIN MRQ CONN AHA CLR CD MTLNK	2017003-001	Cable, ECG, GE, 3/5-Lead, AHA, 3.6 m/12'	Cable, ECG	Reusable	Each	1	Exact	—	
Getinge USA Sales LLC	0012-00-1155-01	CABLE EKG 5 LD SHLD SFTY PT TRNK IABP	2090600-021	Cable, ECG, Mindray Datascope, 5-Lead, AHA, 0.23 m/9"	Cable, ECG	Reusable	Each	1	Exact	—	
AMC	LW-3090R48/5A	SET LDWR 5 LD RDLCNT STD PNCH CLIP UNSHIELD YOKEMATE LWS 48	403751-009	Leadwire, ECG, Reusable, 5-Lead, Radiolucent DIN, AHA	Lead Set, ECG	Reusable	Each	1	Exact	—	
Curbell Medical Products Inc.	MCB-03NB-10AS-0010	SET LEAD DEVICE MED ELECTRICAL-ONLY 3-LEAD ECG CBL FOR	2021141-001	Leadwire Cable, ECG, 3-Lead, Grabbers, AHA	Cable Lead, ECG	Reusable	Each	1	Exact	—	
Curbell Medical Products Inc.	MCB-05NB-10AS-0000	CABLE EKG 5 LD 1 PC MLD PNCH AHA GE 10FT NS LF REUSE	2017003-001 412681-001	Cable, ECG, GE, 3/5-Lead, AHA, 3.6 m/12' Leadwire, ECG, Reusable, 5-Lead, AHA	Cable, ECG Lead Set, ECG	Reusable	Each	1	Clinically Equivalent	We do not offer 5-lead cables/lead combos. Must buy separate cable and lead	
Curbell Medical Products Inc.	RLW-03AW-29AS-0000	LDW REUSEECG LW MULTILINK3LDGRPSNAP29 AHA	411203-001	Leadwire, ECG, Reusable, 3-Lead, AHA, Snaps	Lead Set, ECG	Reusable	Each	1	Exact	—	
Curbell Medical Products Inc.	RLW-03AW-51AS-0000	LEADWIRE EKG 3 LD SNP AHA CLR CD M-LNK 51IN LF	411203-002	Leadwire, ECG, Reusable, 3-Lead, AHA, Snaps	Lead Set, ECG	Reusable	Each	1	Exact	—	
Curbell Medical Products Inc.	RLW-03BK-39AS-0000	LEADWIRE EKG 3 LD PNCH AHA INTELLIVU 39IN RLW	412682-001	Leadwire, ECG, Reusable, 3-Lead, AHA, Grabbers	Lead Set, ECG	Reusable	Each	1	Clinically Equivalent	Shorter length	
Curbell Medical Products Inc.	RLW-03BW-29AS-0000	SET LDWR 3 LD SNP AHA CLR CD M-LNK 29IN NS LF REUSE	412682-001	Leadwire, ECG, Reusable, 3-Lead, AHA, Grabbers	Lead Set, ECG	Reusable	Each	1	Exact	—	

Current Product			AirLife Product							Cross Reference	
Manufacturer	Item	Description	Item	Description	Product Category	Reusable/ Disposable	UOM	QTY	Cross Type	Comments	
Curbell Medical Products Inc.	RLW-03BW-51AS-0000	LDW REUSEECG LW MULTILINK3LDGRPPINCH51 AHA	412682-002	Leadwire, ECG, Reusable, 3-Lead, AHA, Grabbers	Lead Set, ECG	Reusable	Each	1	Exact	—	
Curbell Medical Products Inc.	RLW-05AD-29AS-0000	LDW REUSEECG LW MULTILINK5LEADSNAP29 AHA	412681-001	Leadwire, ECG, Reusable, 5-Lead, AHA, Grabbers	Lead Set, ECG	Reusable	Each	1	Exact	—	
Curbell Medical Products Inc.	RLW-05AW-29AS-0000	LDW REUSEECG LW MULTILINK5LEADGRPSNAP29 AHA	411202-001	Leadwire, ECG, Reusable, 5-Lead, AHA, Snaps	Lead Set, ECG	Reusable	Each	1	Exact	—	
Curbell Medical Products Inc.	RLW-05AW-51AS-0000	LDW REUSEECG LWMULTI-LINK5LEADMLDSNAP51 AHA	411202-002	Leadwire, ECG, Reusable, 5-Lead, AHA, Snaps	Lead Set, ECG	Reusable	Each	1	Exact	—	
Curbell Medical Products Inc.	RLW-05BD-29AS-0000	SET LDWR 5 LD PNCH AHA CLR CD M-LNK 29IN NS LF REUSE	412681-001	Leadwire, ECG, Reusable, 5-Lead, AHA, Grabbers	Lead Set, ECG	Reusable	Each	1	Exact	—	
Curbell Medical Products Inc.	RLW-05BD-51AS-0000	SET LDWR 5 LD PNCH AHA CLR CD M-LNK 51IN NS LF REUSE	412681-002	Leadwire, ECG, Reusable, 5-Lead, AHA, Grabbers	Lead Set, ECG	Reusable	Each	1	Exact	—	
Curbell Medical Products Inc.	RLW-05BM-29AS-0000	SET LDWR 5 LD PNCH AHA CLR CD M-LNK 29IN NS LF REUSE	416467-003	Leadwire, ECG, Reusable, V-Lead, AHA, Grabbers	Lead Set, ECG	Reusable	Each	1	Exact	—	
Curbell Medical Products Inc.	RLW-05BW-29AS-0000	SET LDWR 5 LD PNCH AHA CLR CD M-LNK 29IN NS LF REUSE	412681-001	Leadwire, ECG, Reusable, 5-Lead, AHA, Grabbers	Lead Set, ECG	Reusable	Each	1	Exact	—	
Curbell Medical Products Inc.	RLW-05BW-51AS-0000	LDW REUSEECG LW MULTILINK5LEADGRPPINCH51 AHA	412681-002	Leadwire, ECG, Reusable, 5-Lead, AHA, Grabbers	Lead Set, ECG	Reusable	Each	1	Exact	—	
Curbell Medical Products Inc.	RLW-05BW-MXAS-0000	LDW REUSEECG LW MULTILINK5LDGRPPINCH29-51 MXDA	412681-005	Leadwire, ECG, Reusable, 5-Lead, AHA, Grabbers	Lead Set, ECG	Reusable	Each	1	Exact	—	
Curbell Medical Products Inc.	RLW-06AD-29AS-0000	SET LDWR 6 LD SNP AHA CLR CD M-LNK 29IN NS LF REUSE	421930-001	Leadwire, ECG, Reusable, 6-Lead, AHA, Snaps	Lead Set, ECG	Reusable	Each	1	Exact	—	
Curbell Medical Products Inc.	RLW-06BD-29AS-0000	SET LDWR 6 LD PNCH AHA CLR CD M-LNK 29IN NS LF REUSE	421932-001	Leadwire, ECG, Reusable, 6-Lead, AHA, Grabbers	Lead Set, ECG	Reusable	Each	1	Exact	—	
GE Med Systems Info Tech Equipment	394111-011	SET LDWR 5 LD MTLNK FH TELE	394111-502	Leadwire, ApexPro FH, 5-Lead, Grabber, AHA	Lead Set, ECG	Reusable	Each	1	Exact	—	
Intermed Resources TN LLC	RLW-05BR-51AS-0000	GE 51 IN V LEADS SET GRABBER	416467-001	Leadwire, ECG, Reusable, V-Lead, AHA, Grabbers	Lead Set, ECG	Reusable	Each	1	Exact	—	
Schiller America Inc.	2.410308	SET LEADWIRE ECG WIRE NON-INVASIVE CABLE TRND CR ELECT CV LF	403751-023	Leadwire, ECG, Reusable, 10-Lead, Radiolucent, AHA, Grabbers	Lead Set, ECG	Reusable	Each	1	Exact	—	
Cardinal Health 200 LLC	33103	CABLE ECG KENDALL LD NONSTRL DISP DL AND WIRE SYS 3-LEAD PT	2052133-005	Leadwire, 3-Lead, Grabbers, AHA, GE ApexPro CH	Lead Set, ECG	Disposable	Case	10	Clinically Equivalent	Shorter length	
Cardinal Health 200 LLC	33105	CABLE ECG KENDALL LD NONSTRL DISP DL AND WIRE SYS 5-LEAD PT	2052133-007	Leadwire, 5-Lead, Grabbers, AHA, GE ApexPro CH, 75 cm/29"	Lead Set, ECG	Disposable	Case	10	Clinically Equivalent	Shorter length	
Cardinal Health 200 LLC	33111	LEAD SWITCHING KENDALL ECG NONSTRL REUSABLE DL CABLE	2052133-008	Leadwire, V-Lead, Grabbers, AHA, GE ApexPro CH	Lead Set, ECG	Disposable	Case	10	Clinically Equivalent	Shorter length	

Current Product			AirLife Product							Cross Reference	
Manufacturer	Item	Description	Item	Description	Product Category	Reusable/ Disposable	UOM	QTY	Cross Type	Comments	
Cardinal Health 200 LLC	33114	LEAD EKG 5 LD DRCT CNCT KDL DL LF DISP PHLPS INTELIVU TELEM	2052133-007	Leadwire, 5-Lead, Grabbers, AHA, GE ApexPro CH, 75 cm/29"	Lead Set, ECG	Disposable	Case	10	Clinically Equivalent	Original is Philips Direct Connect, shifted to universal design	
Cardinal Health 200 LLC	33112-	CABLE TELEM GE 5 LD DRCT CNCT GE SYS LF DISP	2052133-007	Leadwire, 5-Lead, Grabbers, AHA, GE ApexPro CH, 75 cm/29"	Lead Set, ECG	Disposable	Case	10	Clinically Equivalent	Shorter length	
Cardinal Health 200 LLC	33126-	SET TELEM APEXPRO 6 LD DRCT CNCT FRQ HOPPING DISP	2086367-004	Leadwire, 6-Lead, Grabbers, AHA, GE ApexPro CH, 75 cm/29"	Lead Set, ECG	Disposable	Case	10	Clinically Equivalent	Original is ApexPro FH Direct Connect, shifted to universal design	
Cardinal Health 200 LLC	33135BT	LEAD ECG 12FT SUSTAINABLE TECHNOLOGIES REPRO NONSTR L DISP	2052133-007	Leadwire, 5-Lead, Grabbers, AHA, GE ApexPro CH, 75 cm/29"	Lead Set, ECG	Disposable	Case	10	Clinically Equivalent	Original is a disposable cable and lead together, no disposable cable needed	
Cardinal Health 200 LLC	33135T	ADAPTER ECG KENDALL NON-INVASIVE CABLE ELECT PHYS MEDICINE	2052133-007	Leadwire, 5-Lead, Grabbers, AHA, GE ApexPro CH, 75 cm/29"	Lead Set, ECG	Disposable	Case	10	Clinically Equivalent	Original is a disposable cable and lead together, no disposable cable needed	
Cardinal Health 200 LLC	33136T	SET LEADWIRE 3FT KENDALL NON-INVASIVE CABLE ELECT PHYS DISP	2052133-007	Leadwire, 5-Lead, Grabbers, AHA, GE ApexPro CH, 75 cm/29"	Lead Set, ECG	Disposable	Case	10	Clinically Equivalent	Shorter length	
Cardinal Health 200 LLC	33145T	CABLE EKG 2 CNCT RDL CNT LD LTWT PLBL RBN DISP	2052133-007	Leadwire, 5-Lead, Grabbers, AHA, GE ApexPro CH, 75 cm/29"	Lead Set, ECG	Disposable	Case	10	Clinically Equivalent	Original is a disposable cable and lead together, no disposable cable needed	
ConMed Corp.	01-3130	BACKPAD EKG ELCTRD ADLT 5% KCL PPR PE FOM ACRYL 7.5X6.5IN	2052133-007	Leadwire, 5-Lead, Grabbers, AHA, GE ApexPro CH, 75 cm/29"	Lead Set, ECG	Disposable	Case	10	Clinically Equivalent	Original is a backpad, new is leadwire	
Curbell Medical Products Inc.	20459	SET LDWR 5 LD AHA CLR CD SNP CONN TRLNK 40IN	2052133-007	Leadwire, 5-Lead, Grabbers, AHA, GE ApexPro CH, 75 cm/29"	Lead Set, ECG	Disposable	Case	10	Clinically Equivalent	Original is a TrueLink connection and SNAP, switched to universal and grabber	
Curbell Medical Products Inc.	DLP-05BD-40AM-0100	LEADWIRE EKG 5 LD PNCH AHA CLR CD ST MTLNK 40IN DISP	2052133-007	Leadwire, 5-Lead, Grabbers, AHA, GE ApexPro CH, 75 cm/29"	Lead Set, ECG	Disposable	Case	10	Clinically Equivalent	Shorter length	

For any questions regarding ECG monitoring, please reach out to your local AirLife Sales Representative.