


HIGH-PERFORMANCE HEATED CIRCUIT TIDAL VOLUME GUIDE



ADULT HEATED CIRCUITS:

ITEM	MINIMUM TIDAL VOLUME	WEIGHT
AH102	120 mL	~ 15 kg
AH202		
AH180		
AH280		

INFANT HEATED CIRCUITS:

ITEM	MAXIMUM TIDAL VOLUME	WEIGHT
AH165	120 mL	~ 15 kg
AH265		

Clinical decisions regarding patient care should be made per patient by a clinical professional. Tidal Volume estimated to be 6-8 mL/kg of ideal body weight.¹

Please refer to your facility's protocol or procedure.



References
 1 Kacmarek R. M., Stoller, J. K., Heuer, A. J., and Egan, D. F. Egan's Fundamentals of Respiratory Care (10th ed.). St. Louis, Mo.: Elsevier/Mosby, 2013.