

HEALTHPARTNERS PRODUCT EVALUATION AIRLIFE™ OPEN MASK




AM-1000 AIRLIFE™ OPEN MASK

- Open oxygen mask featuring dependable oxygen delivery performance and patient comfort
- Provides one solution for all your oxygen delivery needs across your patients' continuum of care, instead of requiring multiple devices for changing flow needs
- Unique shape maintains higher oxygen concentrations while reducing the risk of rebreathing carbon dioxide — all while allowing patients to undertake normal activities like drinking, eating, and talking
- Bifurcated straps eliminate pressure sores on the patient's ears
- Enhanced patient safety that accommodates 1 – 15 LPM flow rates



Scan for
Instructions for Use



HEALTHPARTNERS

IF YOU HAVE ANY QUESTIONS, PLEASE CONTACT:

JUSTIN HALL
TERRITORY MANAGER
952-353-3578

ALETTE JOHNSON
CLINICAL SPECIALIST
248-761-2181

INS-20190-C1 Rev. 1