

AIRFLOW™ ADULT RESUSCITATION MASKS

LATEX FREE | SINGLE PATIENT USE | DISPOSABLE



Ergonomic cushion designed for an optimal seal



Inflation tail valve adjusts cuff pressure for optimal seal



AIRFLOW™ MASKS

- Comes packaged with adult AirFlow resus bags
- Premium soft cushioned mask
- Pliable, ultra-thin wall
- Can be inflated to clinical preference with use of syringe
- Soft seal with exceptional seating characteristics

INS-20192-C1 Rev 1