

RADIOLUCENT LEADWIRES WITH MULTI-LINK™ X2 SYSTEM IEC

SINGLE-SOURCE SOLUTION FOR ALL ECG NEEDS

Unlike many other ECG solution providers, we offer not only reusable and single-patient use portfolios, but also radiolucent leadwires. By standardizing to Multi-Link™ X2, you can choose the Multi-Link leadwires that fit your needs throughout your facility. The Multi-Link X2 ECG System is your single-source solution for all ECG needs.

RADIOLUCENT LEADWIRES FOR X-RAY

As part of the standardized solution, you can use standard reusable Multi-Link X2 ECG cables with the Multi-Link radiolucent leadwires. Choosing your ECG solution for X-ray environment is easy:

SELECT YOUR 3.6 M/12 FT CABLE*:

MONITOR INPUT	SKU	MONITOR BRAND	LEADS	TYPE
	2017003-003	GE HealthCare	3 or 5	Cable
	2017005-003		6	Cable
	2017006-003		12	Cable
	2090200-052	Philips	3 or 5	Cable
	2090200-054		6	Cable
	2090200-056		12	Cable
	2090101-091	Dräger	5 or 6	Adapter
	2090200-062	Mindray Traditional	3 or 5	Cable
	2090200-066		12	Cable
	2090200-072	Spacelabs	5	Cable
	2090200-076		12	Cable

ADD YOUR RADIOLUCENT LEADWIRE*:

	SKU	LEADS	LENGTH
	403751-025	5	150 cm
	403751-033	C	
	403751-026	12	

* These product codes are not available in the United States.



KEY FEATURES

- Multi-Link X2 compatible
- Free of phthalates such as DEHP*
- Does not contain BPA*
- Latex Free
- Designed by ECG experts
- Tested to meet/exceed ECG standards
- Compatibility tested with systems

* These products are not manufactured with natural rubber latex and do not contain DEHP, BBP, DBP or DIBP in concentrations above 0.1% w/w.