


AIRLIFE® OPEN ET AND ET+ MASK

LATEX FREE | SINGLE USE | DISPOSABLE

Introducing AirLife® Open Et and Et+, a new open procedural mask for effective and flexible oxygen delivery with EtCO₂ monitoring. The AirLife Open Et+ Mask is a single use device intended for delivering supplemental oxygen and monitoring exhaled carbon dioxide in non-intubated, spontaneously breathing patients (≥ 20 kg).

FEATURES

- Captures accurate capnography signal while delivering supplemental O₂ so you can monitor respiratory status
- Flexible flow-rates from 1–15 LPM and flush
- Et+ is available with a tethered female/reflective pigtail designed to support EtCO₂ monitoring compatibility throughout the continuum of care
- Distributed flow creates soft oxygen curtain
- Oral access for procedures



AIRLIFE™ OPEN ET AND ET+ MASK

ITEM	DESCRIPTION	SIZE	PK
AM-2000	AirLife™ Open Et Oxygen Mask	Medium	25
AM-2100	AirLife™ Open Et Oxygen Mask	Large	25
AM-2200	AirLife™ Open Et Oxygen Mask	Small	25
AM-3000	AirLife™ Open Et+ Oxygen Mask, Reflective	Medium	25
AM-3100	AirLife™ Open Et+ Oxygen Mask, Reflective	Large	25
AM-3200	AirLife™ Open Et+ Oxygen Mask, Reflective	Small	25