

SUPERNO₂VA[®] ET SYSTEM


Reflective Gas Sampling Line is now included with the SuperNO₂VA[®] Et System



Adult Large Mask



Adult Medium Mask



Compatible with
Microstream Technology



2 L Hyperinflation Bag



Supplemental
O₂ Tubing



Gas Sampling Line with
Male Connectors



SalterSTAT[™] Gas Sampling Line
(Compatible with Microstream technology)

KEY PRODUCT FEATURES

- Includes an additional reflective gas sample line to support ease of use through continuum of care from Pre-Op to PACU
- Supports compatibility with Microstream technology through the use of a reflective gas sampling line
- Packaged in system configuration for supply chain and workflow efficiency

ITEM	DESCRIPTION	MASK SIZE
SSNEL-20	SuperNO ₂ VA Et System with Medium Mask, 2 L Hyperinflation System, O ₂ tubing, Male Gas Sample Line and Reflective Gas Sample Line	Adult Medium
SSNEM-20	SuperNO ₂ VA Et System with Large Mask, 2 L Hyperinflation System, O ₂ Tubing, Male Gas Sample Line and Reflective Gas Sample Line	Adult Large